



**MILLENNIUM  
SURGICAL**

**SURGICALINSTRUMENTS.COM**

UNMATCHED SERVICE • UNLIMITED INSTRUMENTS

822 Montgomery Ave #205  
Narberth, PA 19072  
(800) 600-0428

# Sterilization Processing Instructions for Hospital Instrument Case / Tray Systems

IFU-2017.01.10

This instrument case and tray system has been designed to be compatible with the manufacturer's instructions for steam autoclave sterilization.

To maintain sterility, the instrument case / tray system should be enclosed in an FDA approved sterilization wrap, pouch or rigid sterilization container prior to the sterilization process.

## General Instructions

### Steam Autoclave

Follow manufacturer's recommended processing instructions and ensure that these recommendations are consistent with the surgical instrument manufacturer's processing recommendations (Flash cycles are not recommended). The instrument cases and trays should be dry at the end of the manufacturer's recommended drying cycle. "Wet Packs" requiring additional drying time can occur at the end of the recommended cycle if:

- The instrument system exceeds the 16 - 18 lbs. weight limit.
- The instrument system contains plastic instruments or devices that produce additional moisture requiring additional drying time.

This product has been manufactured using stainless steel, anodized aluminum and medical-grade silicone materials. For further information, contact your Product Consultant.

*Note: Case / Tray Systems require an aluminum safe, neutral pH detergent to avoid faded surface colors and deterioration of the anodized surface.*

IFU 1000-956-002  
Revision D  
Effective 09/2016