



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM


UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Femtosecond Laser Instruments




92-S0118TR: \$287
Thorlakson Adjustable Lid Speculum
Temporal model - 11mm open round blades
with flairs. Stainless Steel.



92-S020020: \$225
Slade-Murdoch Eye Speculum
Self-Lock w/Control Screw Blade
conforms to the curvature of the
femtosecond laser. 20mm open round
blades with flairs. Stainless Steel.



92-H0359L: \$245
Chang Combination Chopper
A combination of an altered MicroFinger and
a quick chop tip.



92-H0637: \$163
Nagahara-Ahn Small Pupil Chopper
1.25mm inner cutting edge (for small pupils
or floppy iris). Notch on the blunt side (used
to retract pupil).

Visit our website for a complete line of Femtosecond Laser Instruments
www.surgicalinstruments.com



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428


Femtosecond Laser Instruments



92-F25213: \$375

Tjia-Akahoshi Utrata Capsulorhexis Forceps

Long, thin blades for use in 2mm or smaller incisions.
Vaulted jaws. Contours to the anterior capsule.



92-816172: \$187

Intralase Flap Lifter

Intra Flap Lifter: Slightly-pointed tip, angled, bend to tip:
1.5mm. Retreatment Spatula: Z-shaped tip, 12mm wide,
cylindrical spatula w/ slightly pointed tip, 106mm length.



92-H0279: \$292

Thorlakson Flap Lifter

Combines a modified Sinskey Hook and a Bi-Convex
Spatula with beveled notches, angled / vaulted, 12mm
from bend to tip, round knurled handle, 133mm overall
length. Stainless Steel.



92-H0499: \$222

Koch Spatula and Nagahara Chopper

Overall length: 130mm. Cutting edge: 1.5mm. Cutting
edge is 60 degrees (cuts toward phaco tip). Manipulator
pushes, rotates, lifts, and splits - serves as protective
barrier between phaco tip and posterior capsule.

92-H0706



92-H0706

Femtosecond Double End Spatula.
Spear end is angled 30 degrees.

Visit our website for a complete line of Femtosecond Laser Instruments
www.surgicalinstruments.com



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Femtosecond Laser Instruments

Chang Combination Chopper



92-H0359L

92-H0359L

A combination of an altered MicroFinger and a quick chop tip.

Murdoch Eye Speculum



3-4094

3-4094

Murdoch Eye Speculum

Nagahara-Ahn Small Pupil Chopper



92-H0637

92-H0637

1.25mm inner cutting edge (for small pupils or floppy iris). Notch on the blunt side (used to retract pupil).

Akahoshi Utrata Capsulorhexis Forceps



92-F2521

92-F2521

Long, thin blades for use in 2mm or smaller incisions. 11mm blade length. Sharp tip.

Tjia-Akahoshi Utrata Capsulorhexis Forceps



92-F25213

92-F25213

Long, thin blades for use in 2mm or smaller incisions. Vaulted jaws. Contours to the anterior capsule.

Intralase Flap Lifter



92-H0356

92-H0356

Double ended.

Thorlakson Flap Lifter



92-H0279

92-H0279

Combines a modified Sinsky Hook and a Bi-Convex Spatula with beveled notches, angled / vaulted, 12mm from bend to tip, round knurled handle, 133mm overall length. Stainless Steel.

Koch Spatula and Nagahara Chopper



92-H0499

92-H0499

Overall length: 130mm. Cutting edge: 1.5mm. Cutting edge is 60 degrees (cuts toward phaco tip). Manipulator pushes, rotates, lifts, and splits - serves as protective barrier between phaco tip and posterior capsule.



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Femtosecond Laser Instruments

Chang Combination Chopper

92-H0359L: A combination of an altered MicroFinger and a quick chop tip.

Nagahara-Ahn Small Pupil Chopper

92-H0637
1.25mm inner cutting edge (for small pupils or floppy iris). Notch on the blunt side (used to retract pupil).

Tjia-Akahoshi Utrata Capsulorhexis Forceps

92-F25213
Long, thin blades for use in 2mm or smaller incisions. Vaulted jaws. Contours to the anterior capsule.

Thorlakson Flap Lifter

92-H0279
Combines a modified Sinskey Hook and a Bi-Convex Spatula with beveled notches, angled / vaulted, 12mm from bend to tip, round knurled handle, 133mm overall length. Stainless Steel.

Murdoch Eye Speculum

3-4094: Murdoch Eye Speculum

Akahoshi Utrata Capsulorhexis Forceps

92-F2521
Long, thin blades for use in 2mm or smaller incisions. 11mm blade length. Sharp tip.

Koch Spatula and Nagahara Chopper

92-H0499
Overall length: 130mm. Cutting edge: 1.5mm. Cutting edge is 60 degrees (cuts toward phaco tip). Manipulator pushes, rotates, lifts, and splits - serves as protective barrier between phaco tip and posterior capsule.

Intralase Flap Lifter

92-H0356: Double ended.



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM


UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Femtosecond Laser Instruments



92-H0359L



3-4094




92-H0637



92-F2521



92-F25213




92-H0356



92-H0279



92-H0499



92-H0706A



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Femtosecond Laser Instruments

Product #	Item	Description	Price
92-F25133	Utrata MICS Capsulorhexis Forceps	Jaws set to open 1.8mm from tip to bend. Angled 12mm from bend to tip, extremely thin triangular grasping tips, serrated flat handle, 80mm overall length. Stainless Steel.	
92-S020020	Slade-Murdock Speculum	Self-Lock w/Control Screw Blade conforms to the curvature of the femtosecond laser. 20mm open round blades with flairs. Stainless Steel.	
92-4421H27	Nichamin Hydrodissection Cannula	27ga flattened tip, angled / vaulted, 12mm from bend to tip, 22mm overall length.	
92-H0304	Slade Femtosecond Spatula	1mm duck-bill flat tip, angled / semi-vaulted. 10mm from bend to tip. Round knurled handle. 114mm overall length. Stainless Steel. Used to open the main and side port incisions made by the femtosecond laser.	
92-S0118TR	Thorlakson Adjustable Lid Speculum	Temporal model - 11mm open round blades with flairs. Stainless Steel.	
92-816172	TTI Intra Flap Lifter & Retreatment Spatula	Double-ended. Intra Flap Lifter: Delicate slightly-pointed tip, angled, 1.5mm from bend to tip. Retreatment Spatula: Unique z-shaped tip, 12mm wide cylindrical spatula with slightly pointed tip, round knurled handle, 106mm overall length. Stainless Steel.	
92-C0469	Slade Femtosecond Hydrodissection Cannula	23ga, angled / vaulted, 14mm from bend to tip, vertically flattened tip with a posterior flow, 20mm overall length.	
3-4072	Mellinger Eye Speculum	With Fenestrated Blades	



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Femtosecond Laser Instruments

