



**MILLENNIUM
SURGICAL**


SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS


822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Retina / Vitrectomy

23 and 20 Gauge Vitreo-Retinal Scissors and Forceps




High-precision titanium tips on the new titanium handle for vitreoretinal surgery



By exploring advanced concepts and embracing new technology, Duckworth and Kent has created an innovative range of vitrectomy instruments utilising the new ergonomic handle design. From the surgeon's perspective, the transition to the new style will feel both seamless and rewarding.

The new design features detachable tip heads for cleaning purposes. These tips can be screwed onto the threaded handle, a mechanism which eliminates the need for alignment when engaging the head onto the instrument handle.


The new retinal instruments are supplied complete (handle and head) in own specifically designed sterilising tray (as shown above).




Protective Cleaning Cover




The new cover is used to protect the detachable tip heads during cleaning and reprocess cycles. A luer flushing adaptor at the back of the cover enables thorough cleaning and disinfection of the fine delicate tips.



Tip head screws on



Tip head pivots then locks into place securely



Syringe or other luer device back flushes tip head



MILLENNIUM
SURGICAL

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Retina / Vitrectomy

23 and 20 Gauge Vitreo-Retinal Scissors and Forceps

Scissors




87-1842N



Curved Scissors

NEW tip design 23G

87-1841N



Straight Scissors

23G


87-1841BN



Straight Blunt Scissors

23G

87-18411N



Straight Scissors

20G

Forceps



87-2879N



Asymmetrical Forceps

NEW tip design 23G

87-2877N



End Gripping Forceps
ILM Peeling

23G


87-28781N



Delicate Serrated
Forceps

23G


87-2868N



Smooth Rounded Tip
Forceps

20G


87-2871N



End Gripping Forceps
ILM Peeling

20G

87-2878N



Delicate Serrated
Forceps

20G

87-2885N



Rassam Membrane
Picking Forceps

20G