



**MILLENNIUM
SURGICAL**


SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Endoscope Holder/Positioner

The Goose Neck Retractor System enhances is mounted on the accessory rails of the operating tables by means of the Steril-Clamp™. Flexibility and maneuverability are facilitated by fifteen individual joints, based on the well established Goose-Neck principle. The different types of instrument holders that we have available add many unique features to this retractor system.



72-1509060 (A): Goose neck retractor system with self-retaining sleeve clamp with articulating joint for holding telescopes. Includes 72-1508056 (Steril-Clamp™ for attaching to table).

72-1509061 (B): Goose neck retractor system with universal mounting clamp for holding all types of instruments. Includes 72-1508056 (Steril-Clamp™ for attaching to table).

FEATURES:

- All components made of a high-grade stainless steel.
- The complete system, including the Steril-Clamp™, may be flash autoclaved.
- The self-retaining joints rotate and twist in every desirable position.
- The instrument shaft fits into table accessory clamps such as those used for stirrups and leg holders, as well as all of the other over-the-drape clamps with an opening of 5/8".
- The vertical shaft (5/8" Ø) may be lowered or heightened within the clamp for flexible height adjustment.
- Two different models with different instrument holders are available.
- They are economically priced.



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Endoscope Holder/Positioner

Scope Holding Device

Video endoscopic surgery has become a standard procedure in every hospital & surgery center throughout the country and in most parts of the world. The Endo-Fix™ Scope Holding Device is based on over 30 years of experience in the design, manufacturing and marketing of the well known Retract-Robot™ and the Endoscopic Fixation Device.



**Simple articulation and precise holding
with spring loaded pressure release joints:
The most economical-priced system made in the
U.S.A.**

72-1509060: Goose neck retractor system with self-retaining sleeve clamp with articulating joint. Includes 72-1508056 (Steril-Clamp™ for attaching to table.

FEATURES:

- Self-retaining and positioning.
- Will hold laparoscope with camera in any desirable position.
- Simplified operation and reposition without opening/tightening the friction joints.
- Will hold other surgical instruments, such as retractors, graspers etc.
- Unique alligator-type instrument holder for endoscopes and instruments.
- The instruments may be exchanged without violating sterile requirements.
- Mounts on the accessory bar of any operating room table by means of the autoclavable steril-clamp™.
- Easy autoclavable without disassembly. Just fold, wrap and place into conventional autoclave.
- Economically priced, but optimally designed for performance and durability.
- All parts are made out of medical grade stainless steel and are completely manufactured in the U.S.A.