



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Laparoscopic Instruments

Reusable 32 and 45 cm length laparoscopic instruments

- 32 and 45 cm lengths accommodate standard and advanced laparoscopic procedures, including Gynecology & Bariatrics.
- Rotating shafts with 360 degree rotation, one finger controls the tip placement.
- A wide variety of handle styles are available to accommodate a wide variety of surgeon preference.
- Extensive selection of tip styles.
- Offered in two styles to meet your cleaning requirements.
- Modular Laparoscopic instruments can be quickly disassembled for cleaning and visual inspection. Processing technicians spend less time removing bio-burden and the risk of HAI decreases when they can be certain that bio-burden is not trapped in the instrument. Easily reassembled for the surgeon prior to use in surgery.

Single Piece with Flush-Port for quick flushing of bio-burden



Features:

- Durability and Longevity - Stainless steel construction, medical grade polymer coating guarantee long life, easy maintenance and safety.
- Precise & Ergonomic Handling - Jaw patterns are designed and machined for precision grasping, dissecting or cutting. Large-finger handles with tang and large diameter rotating drum with stainless steel roller bearings give you precise control and comfort while reducing hand fatigue.
- Optimal Cleaning & Maintenance - Flush ports and expanded push-rod clearance on non-modular instruments, minimum components and solid connections on modular instruments assures easy maintenance.
- Excellent Return on Investment - Completely reusable and fully-autoclavable.

Handle Styles



TYPE 1:
INSULATED WITH RF-
POST



TYPE 2:
INSULATED W/ RF-
POST & RATCHET



TYPE 3:
NON-INSULATED
STAINLESS STEEL



TYPE 4:
NON-INSULATED
STAINLESS STEEL W/
RATCHET



TYPE 5:
LIGHTWEIGHT
CARBON-FIBER W/ RF
POST



TYPE 6:
LIGHTWEIGHT
CARBON-FIBER
w/ RF-POST &
RATCHET



TYPE 7:
LIGHTWEIGHT
CARBON-FIBER
W/ RF POST,
RATCHET /
RATCHET
DEFEATING



TYPE 8:
LIGHTWEIGHT CARBON
-FIBER W/ RF POST,
FRENCH ON/OFF
RATCHET



TYPE 9:
LIGHTWEIGHT
CARBON-FIBER W/ RF
POST, W/ LARGE
FINGER LOOPS



TYPE 10:
LIGHTWEIGHT
CARBON-FIBER W/ RF
POST, RATCHET /
RATCHET DEFEATING
W/ LARGE FINGER
LOOPS

Ordering Example: 72-5225501C1

TIP

STYLE HANDLE



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

5 MM Ø Dissectors & Graspers

		Rigid One-Piece with Flush Port	Rotatable One-Piece with Flush Port	
	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>	
	5MM Ø, MICRO REDDICK-OLSEN GRASPER/DISSECTOR	33 CM	52255-01A	52255-01B
		36 CM	52236-01A	52236-01B
		45 CM	52155-01A	52155-01B
		24 CM	52224-01A	52224-01B
	5MM Ø, STANDARD GRASPER/DISSECTOR WITH SPOON	33 CM	52255-02A	52255-02B
		36 CM	52236-02A	52236-02B
		45 CM	52155-02A	52155-02B
		24 CM	52224-02A	52224-02B
	5MM Ø, STANDARD GRASPER/DISSECTOR W/ SPOON, LONG 20 MM JAW	52255-01A	52255-52A	52255-52B
		52236-01A	52236-52A	52236-52B
		52155-01A	52155-52A	52155-52B
		52221-01A	52224-52A	52224-52B
	5MM Ø, STANDARD GRASPER/DISSECTOR	33 CM	52255-03A	52255-03B
		36 CM	52236-03A	52236-03B
		45 CM	52155-03A	52155-03B
		24 CM	52224-03A	52224-03B
	5MM Ø, BULLET-NOSE GRASPER/DISSECTOR WITH SPOON	33 CM	52255-04A	52255-04B
		36 CM	52236-04A	52236-04B
		45 CM	52155-04A	52155-04B
		24 CM	52224-04A	52224-04B
	5MM Ø, TAPERED DOLPHIN-NOSE GRASPER/DISSECTOR	33 CM	52255-05A	52255-05B
		36 CM	52236-05A	52236-05B
		45 CM	52155-05A	52155-05B
		24 CM	52224-05A	52224-05B
	5MM Ø, TAPERED DOLPHIN-NOSE GRASPER/DISSECTOR W/ SPOON	33 CM	52255-105A	52255-105B
		36 CM	52236-105A	52236-105B
		45 CM	52155-105A	52155-105B
		24 CM	52224-105A	52224-105B
	5MM Ø, MARYLAND GRASPER/DISSECTOR, 18MM LENGTH TIP	33 CM	52255-06A	52255-06B
		36 CM	52236-06A	52236-06B
		45 CM	52155-06A	52155-06B
		24 CM	52224-06A	52224-06B
	5MM Ø, TAPERED MARYLAND GRASPER/DISSECTOR, 18MM LENGTH TIP	33 CM	52255-99A	52255-99B
		36 CM	52236-99A	52236-99B
		45 CM	52155-99A	52155-99B
		24 CM	52224-99A	52224-99B
	5MM Ø, SHORT MARYLAND GRASPER/DISSECTOR, 15MM LENGTH TIP	33 CM	52255-91A	52255-91B
		36 CM	52236-91A	52236-91B
		45 CM	52155-91A	52155-91B
		24 CM	52224-91A	52224-91B



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

5 MM Ø Dissectors & Graspers

Rigid One-Piece with Flush Port Rotatable One-Piece with Flush Port

	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
	33 CM	52255-07A	52255-07B
	36 CM	52236-07A	52236-07B
	45 CM	52155-07A	52155-07B
	24 CM	52224-07A	52224-07B
	33 CM	52255-59A	52255-59B
	36 CM	52236-59A	52236-59B
	45 CM	52155-59A	52155-59B
	24 CM	52224-59A	52224-59B
	33 CM	52255-08A	52255-08B
	36 CM	52236-08A	52236-08B
	45 CM	52155-08A	52155-08B
	24 CM	52224-08A	52224-08B
	33 CM	52255-09A	52255-09B
	36 CM	52236-09A	52236-09B
	45 CM	52155-09A	52155-09B
	24 CM	52224-09A	52224-09B
	33 CM	52255-10A	52255-10B
	36 CM	52236-10A	52236-10B
	45 CM	52155-10A	52155-10B
	24 CM	52224-10A	52224-10B
	33 CM	52255-11A	52255-11B
	36 CM	52236-11A	52236-11B
	45 CM	52155-11A	52155-11B
	24 CM	52224-11A	52224-11B
	33 CM	52255-53A	52255-53B
	36 CM	52236-53A	52236-53B
	45 CM	52155-53A	52155-53B
	24 CM	52224-53A	52224-53B
	33 CM	52255-12A	52255-12B
	36 CM	52236-12A	52236-12B
	45 CM	52155-12A	52155-12B
	24 CM	52224-12A	52224-12B
	33 CM	52255-13A	52255-13B
	36 CM	52236-13A	52236-13B
	45 CM	52155-13A	52155-13B
	24 CM	52224-13A	52224-13B
	33 CM	52255-14A	52255-14B
	36 CM	52236-14A	52236-14B
	45 CM	52155-14A	52155-14B
	24 CM	52224-14A	52224-14B



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

5 MM Ø Dissectors & Graspers

		Rigid One-Piece with Flush Port	Rotatable One-Piece with Flush Port
	<i>Length</i>	<i>Style A:</i>	<i>Style B:</i>
	33 CM	52255-15A	52255-15B
	36 CM	52236-15A	52236-15B
	45 CM	52155-15A	52155-15B
	24 CM	52224-15A	52224-15B
5MM Ø, BABCOCK GRASPER, FENESTRATED TIP	33 CM	52255-16A	52255-16B
	36 CM	52236-16A	52236-16B
	45 CM	52155-16A	52155-16B
	24 CM	52224-16A	52224-16B
5MM Ø, BABCOCK GRASPER, AXIS STYLE, FENESTRATED TIP	33 CM	52255-116A	52255-116B
	36 CM	52236-116A	52236-116B
	45 CM	52155-116A	52155-116B
	24 CM	52224-116A	52224-116B
5MM Ø, BABCOCK GRASPER- SHORT: DELICATE APPLICATION	33 CM	52255-17A	52255-17B
	36 CM	52236-17A	52236-17B
	45 CM	52155-17A	52155-17B
	24 CM	52224-17A	52224-17B
5MM Ø, ALLIS GRASPER, SHORT VERSION, SINGLE ACTION	33 CM	52255-18A	52255-18B
	36 CM	52236-18A	52236-18B
	45 CM	52155-18A	52155-18B
	24 CM	52224-18A	52224-18B
5MM Ø, ALLIS GRASPER, SHORT VERSION, 3 X 4 TEETH, SERRATED, DOUBLE- ACTION	33 CM	52255-19A	52255-19B
	36 CM	52236-19A	52236-19B
	45 CM	52155-19A	52155-19B
	24 CM	52224-19A	52224-19B
5MM Ø, ATRAUMATIC ALLIS GRASPER, SHORT VERSION, SERRATED W/ DOUBLE ROW OF TEETH, DOUBLE-ACTION	33 CM	52255-119A	52255-119B
	36 CM	52236-119A	52236-119B
	45 CM	52155-119A	52155-119B
	24 CM	52224-119A	52224-119B
5MM Ø, ATRAUMATIC ALLIS GRASPER, LONG VERSION, SERRATED W/ DOUBLE ROW OF TEETH, DOUBLE-ACTION	33 CM	52255-219A	52255-219B
	36 CM	52236-219A	52236-219B
	45 CM	52155-219A	52155-219B
	24 CM	52224-219A	52224-219B
5MM Ø, ALLIS GRASPER, MEDIUM VERSION, (WECK STYLE)	33 CM	52255-20A	52255-20B
	36 CM	52236-20A	52236-20B
	45 CM	52155-20A	52155-20B
	24 CM	52224-20A	52224-20B
5MM Ø, ALLIS GRASPER, STANDARD VERSION	33 CM	52255-21A	52255-21B
	36 CM	52236-21A	52236-21B
	45 CM	52155-21A	52155-21B
	24 CM	52224-21A	52224-21B



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

5 MM Ø Dissectors & Graspers

Rigid
One-Piece
with Flush Port

Rotatable
One-Piece
with Flush Port

	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
	33 CM	52255-57A	52255-57B
	36 CM	52236-57A	52236-57B
	45 CM	52155-57A	52155-57B
	24 CM	52224-57A	52224-57B
5MM Ø, DEBAKEY GRASPER, 20 MM JAW LENGTH	33 CM	52255-22A	52255-22B
	36 CM	52236-22A	52236-22B
	45 CM	52155-22A	52155-22B
	24 CM	52224-22A	52224-22B
5MM Ø, DEBAKEY GRASPER, 40 MM JAW LENGTH	33 CM	52255-23A	52255-23B
	36 CM	52236-23A	52236-23B
	45 CM	52155-23A	52155-23B
	24 CM	52224-23A	52224-23B
5MM Ø, HUNTER/ FENESTRATED DEBAKEY GRASPER, 20 MM JAW LENGTH	33 CM	52255-58A	52255-58B
	36 CM	52236-58A	52236-58B
	45 CM	52155-58A	52155-58B
	24 CM	52224-58A	52224-58B
5MM Ø, HUNTER/ FENESTRATED DEBAKEY GRASPER, 40 MM JAW LENGTH	33 CM	52255-159A	52255-159B
	36 CM	52236-159A	52236-159B
	45 CM	52155-159A	52155-159B
	24 CM	52224-159A	52224-159B
5MM Ø, HUNTER/ FENESTRATED DEBAKEY GRASPER, 60 MM JAW LENGTH	33 CM	52255-160A	52255-160B
	36 CM	52236-160A	52236-160B
	45 CM	52155-160A	52155-160B
	24 CM	52224-160A	52224-160B
5MM Ø, DEBAKEY CURVED GRASPER, w/ ATRAUMATIC TEETH, DOUBLE-ACTION	33 CM	52255-46A	52255-46B
	36 CM	52236-46A	52236-46B
	45 CM	52155-46A	52155-46B
	24 CM	52224-46A	52224-46B
5MM Ø, GLASSMAN ATRAUMATIC GRASPER	33 CM	52255-24A	52255-24B
	36 CM	52236-24A	52236-24B
	45 CM	52155-24A	52155-24B
	24 CM	52224-24A	52224-24B
5MM Ø, ATRAUMATIC GRASPER, CROSS SERRATED	33 CM	52255-25A	52255-25B
	36 CM	52236-25A	52236-25B
	45 CM	52155-25A	52155-25B
	24 CM	52224-25A	52224-25B
5MM Ø, DOUBLE CUP FUNDUS MAXI GRASPER	33 CM	52255-26A	52255-26B
	36 CM	52236-26A	52236-26B
	45 CM	52155-26A	52155-26B
	24 CM	52224-26A	52224-26B



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

5 MM Ø Dissectors & Graspers



		Rigid One-Piece with Flush Port	Rotatable One-Piece with Flush Port
	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
	33 CM	52255-326A	52255-326B
	36 CM	52236-326A	52236-326B
	45 CM	52155-326A	52155-326B
	24 CM	52224-326A	52224-326B
5MM Ø, "MCKERNAN" STYLE FUNDUS MAXI GRASPER	33 CM	52255-126A	52255-126B
	36 CM	52236-126A	52236-126B
	45 CM	52155-126A	52155-126B
	24 CM	52224-126A	52224-126B
5MM Ø, DOUBLE CUP FENESTRATED FUNDUS MAXI GRASPER	33 CM	52255-226A	52255-226B
	36 CM	52236-226A	52236-226B
	45 CM	52155-226A	52155-226B
	24 CM	52224-226A	52224-226B
5MM Ø, ENDO CLINCH FENESTRATED GRASPER	33 CM	52255-27A	52255-27B
	36 CM	52236-27A	52236-27B
	45 CM	52155-27A	52155-27B
	24 CM	52224-27A	52224-27B
5MM Ø, RETRO FACING RAT-TOOTH GRASPER w/ ATRAUMATIC CUP	33 CM	52255-28A	52255-28B
	36 CM	52236-28A	52236-28B
	45 CM	52155-28A	52155-28B
	24 CM	52224-28A	52224-28B
5MM Ø, ALLIGATOR GRASPER, 10MM JAW, RETRO FACING TEETH, DBL-ACTION	33 CM	52255-29A	52255-29B
	36 CM	52236-29A	52236-29B
	45 CM	52155-29A	52155-29B
	24 CM	52224-29A	52224-29B
5MM Ø, ALLIGATOR GRASPER, 20MM JAW , RETRO FACING TEETH, DBL-ACTION	33 CM	52255-129A	52255-129B
	36 CM	52236-129A	52236-129B
	45 CM	52155-129A	52155-129B
	24 CM	52224-129A	52224-129B
5MM Ø, ALLIGATOR GRASPER, RETRO FACING TEETH, SINGLE-ACTION	33 CM	52255-30A	52255-30B
	36 CM	52236-30A	52236-30B
	45 CM	52155-30A	52155-30B
	24 CM	52224-30A	52224-30B
5MM Ø, CUP SHAPED GRASPER WITH RETRO FACING TEETH, DBL-ACTION	33 CM	52255-58A	52255-58B
	36 CM	52236-58A	52236-58B
	45 CM	52155-58A	52155-58B
	24 CM	52224-58A	52224-58B
5MM Ø, WAVE GRASPER, LONG BEVELED EDGE JAW, DBL-ACTION	33 CM	52255-32A	52255-32B
	36 CM	52236-32A	52236-32B
	45 CM	52155-32A	52155-32B
	24 CM	52224-32A	52224-32B



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

5 MM Ø Dissectors & Graspers



Rigid
One-Piece
with Flush Port

Rotatable
One-Piece
with Flush Port

	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
	33 CM	52255-33A	52255-33B
	36 CM	52236-33A	52236-33B
	45 CM	52155-33A	52155-33B
	24 CM	52224-33A	52224-33B
	33 CM	52255-36A	52255-36B
	36 CM	52236-36A	52236-36B
	45 CM	52155-36A	52155-36B
	24 CM	52224-36A	52224-36B
	33 CM	52255-136A	52255-136B
	36 CM	52236-136A	52236-136B
	45 CM	52155-136A	52155-136B
	24 CM	52224-136A	52224-136B
	33 CM	52255-34A	52255-34B
	36 CM	52236-34A	52236-34B
	45 CM	52155-34A	52155-34B
	24 CM	52224-34A	52224-34B
	33 CM	52255-35A	52255-35B
	36 CM	52236-35A	52236-35B
	45 CM	52155-35A	52155-35B
	24 CM	52224-35A	52224-35B
	33 CM	52255-134A	52255-134B
	36 CM	52236-134A	52236-134B
	45 CM	52155-134A	52155-134B
	24 CM	52224-134A	52224-134B
	33 CM	52255-135A	52255-135B
	36 CM	52236-135A	52236-135B
	45 CM	52155-135A	52155-135B
	24 CM	52224-135A	52224-135B
	33 CM	52255-37A	52255-37B
	36 CM	52236-37A	52236-37B
	45 CM	52155-37A	52155-37B
	24 CM	52224-37A	52224-37B
	33 CM	52255-38A	52255-38B
	36 CM	52236-38A	52236-38B
	45 CM	52155-38A	52155-38B
	24 CM	52224-38A	52224-38B
	33 CM	52255-39A	52255-39B
	36 CM	52236-39A	52236-39B
	45 CM	52155-39A	52155-39B
	24 CM	52224-39A	52224-39B



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

5 MM Ø Dissectors & Graspers

	<u>Length</u>	Rigid	Rotatable
		One-Piece with Flush Port	One-Piece with Flush Port
		<u>Style A:</u>	<u>Style B:</u>
	33 CM	52255-40A	52255-40B
	36 CM	52236-40A	52236-40B
	45 CM	52155-40A	52155-40B
	24 CM	52224-40A	52224-40B
5MM Ø, JOHAN FENESTRATED GRASPER, 22 MM SINGLE-ACTION FLAT JAWS	33 CM	52255-41A	52255-41B
	36 CM	52236-41A	52236-41B
	45 CM	52155-41A	52155-41B
	24 CM	52224-41A	52224-41B
5MM Ø, JOHAN FENESTRATED GRASPER, 45 MM SGL-ACTION FLAT JAWS	33 CM	52255-42A	52255-42B
	36 CM	52236-42A	52236-42B
	45 CM	52155-42A	52155-42B
	24 CM	52224-42A	52224-42B
5MM Ø, KOCHER GRASPER, 1 X 2 TEETH, SINGLE-ACTION FLAT JAWS	33 CM	52255-43A	52255-43B
	36 CM	52236-43A	52236-43B
	45 CM	52155-43A	52155-43B
	24 CM	52224-43A	52224-43B
5MM Ø, KOCHER GRASPER, 1 X 2 TEETH, DOUBLE-ACTION (WECK TYPE)	33 CM	52255-44A	52255-44B
	36 CM	52236-44A	52236-44B
	45 CM	52155-44A	52155-44B
	24 CM	52224-44A	52224-44B
5MM Ø, SCHROEDER TENACULUM GRASPER, SINGLE TEETH	33 CM	52255-45A	52255-45B
	36 CM	52236-45A	52236-45B
	45 CM	52155-45A	52155-45B
	24 CM	52224-45A	52224-45B
5MM Ø, SCHROEDER TENACULUM GRASPER, SINGLE TEETH w/ EXTRA GUIDE TOOTH	33 CM	52255-94A	52255-94B
	36 CM	52236-94A	52236-94B
	45 CM	52155-94A	52155-94B
	24 CM	52224-94A	52224-94B
5MM Ø, BLUNT ATRAUMATIC GRASPING FORCEPS, DOUBLE-ACTION	33 CM	52255-47A	52255-47B
	36 CM	52236-47A	52236-47B
	45 CM	52155-47A	52155-47B
	24 CM	52224-47A	52224-47B
5MM Ø, BLUNT ATRAUMATIC GRASPER W/ CROSS-SERRATED TIPS, DBL-ACTION	33 CM	52255-56A	52255-56B
	36 CM	52236-56A	52236-56B
	45 CM	52155-56A	52155-56B
	24 CM	52224-56A	52224-56B
5MM Ø, BULLET NOSE GRASPER/DISSECTOR, TC, MICRO DIAMOND JAW PATTERN DOUBLE-ACTION	33 CM	52255-48A	52255-48B
	36 CM	52236-48A	52236-48B
	45 CM	52155-48A	52155-48B
	24 CM	52224-48A	52224-48B



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

5 MM Ø Dissectors & Graspers

		Rigid One-Piece with Flush Port	Rotatable One-Piece with Flush Port
	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
	33 CM	52255-49A	52255-49B
	36 CM	52236-49A	52236-49B
	45 CM	52155-49A	52155-49B
	24 CM	52224-49A	52224-49B
	33 CM	52255-51A	52255-51B
	36 CM	52236-51A	52236-51B
	45 CM	52155-51A	52155-51B
	24 CM	52224-51A	52224-51B
	33 CM	52255-54A	52255-54B
	36 CM	52236-54A	52236-54B
	45 CM	52155-54A	52155-54B
	24 CM	52224-54A	52224-54B
	33 CM	52255-55A	52255-55B
	36 CM	52236-55A	52236-55B
	45 CM	52155-55A	52155-55B
	24 CM	52224-55A	52224-55B
	33 CM	52255-90A	52255-90B
	36 CM	52236-90A	52236-90B
	45 CM	52155-90A	52155-90B
	24 CM	52224-90A	52224-90B
	33 CM	52255-91A	52255-91B
	36 CM	52236-91A	52236-91B
	45 CM	52155-91A	52155-91B
	24 CM	52224-91A	52224-91B
	33 CM	52255-144A	52255-144B
	36 CM	52236-144A	52236-144B
	45 CM	52155-144A	52155-144B
	24 CM	52224-144A	52224-144B
	33 CM	52255-145A	52255-145B
	36 CM	52236-145A	52236-145B
	45 CM	52155-145A	52155-145B
	24 CM	52224-145A	52224-145B



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

5MM Biopsy Forceps & Scissors

		Rigid One-Piece with Flush Port	Rotatable One-Piece with Flush Port
	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
	33 CM	52255-60A	52255-60B
	36 CM	52236-60A	52236-60B
	45 CM	52155-60A	52155-60B
	24 CM	52224-60A	52224-60B
5MM Ø, MICRO TIP SCISSORS, STRAIGHT, SINGLE-ACTION	33 CM	52255-61A	52255-61B
	36 CM	52236-61A	52236-61B
	45 CM	52155-61A	52155-61B
	24 CM	52224-61A	52224-61B
5MM Ø, MICRO TIP SCISSORS, CURVED, SINGLE-ACTION	33 CM	52255-62A	52255-62B
	36 CM	52236-62A	52236-62B
	45 CM	52155-62A	52155-62B
	24 CM	52224-62A	52224-62B
5MM Ø, METZENBAUM SCISSORS, STRAIGHT, LONG BLADE (20MM TIP)	33 CM	52255-63A	52255-63B
	36 CM	52236-63A	52236-63B
	45 CM	52155-63A	52155-63B
	24 CM	52224-63A	52224-62B
5MM Ø, METZENBAUM SCISSORS, STRAIGHT, SHORT BLADE (12MM TIP)	33 CM	52255-64A	52255-64B
	36 CM	52236-64A	52236-64B
	45 CM	52155-64A	52155-64B
	24 CM	52224-64A	52224-64B
5MM Ø, METZENBAUM SCISSORS, CURVED, LONG BLADE (18MM TIP)	33 CM	52255-65A	52255-65B
	36 CM	52236-65A	52236-65B
	45 CM	52155-65A	52155-65B
	24 CM	52224-65A	52224-65B
5MM Ø, METZENBAUM SCISSORS, HARDENED STEEL, CURVED BLADE (16MM TIP)	33 CM	52255-66A	52255-66B
	36 CM	52236-66A	52236-66B
	45 CM	52155-66A	52155-66B
	24 CM	52224-66A	52224-66B
5MM Ø, METZENBAUM SCISSORS, CURVED, SHORT BLADE (11MM TIP)	33 CM	52255-166A	52255-166B
	36 CM	52236-166A	52236-166B
	45 CM	52155-166A	52155-166B
	24 CM	52224-166A	52224-166B
5MM Ø, PERITONEAL SCISSORS, SINGLE-ACTION	33 CM	52255-67A	52255-67B
	36 CM	52236-67A	52236-67B
	45 CM	52155-67A	52155-67B
	24 CM	52224-67A	52224-67B
5MM Ø, HOOK SCISSORS, TS, EXTRA HARD STEEL	33 CM	52255-76A	52255-76B
	36 CM	52236-76A	52236-76B
	45 CM	52155-76A	52155-76B
	24 CM	52224-76A	52224-76B



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

5MM Biopsy Forceps & Scissors

		Rigid One-Piece with Flush Port	Rotatable One-Piece with Flush Port
	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
	33 CM	52255-77A	52255-77B
	36 CM	52236-77A	52236-77B
	45 CM	52155-77A	52155-77B
	24 CM	52224-77A	52224-77B
	33 CM	52255-78A	52255-78B
	36 CM	52236-78A	52236-78B
	45 CM	52155-78A	52155-78B
	24 CM	52224-78A	52224-78B
	33 CM	52255-79A	52255-79B
	36 CM	52236-79A	52236-79B
	45 CM	52155-79A	52155-79B
	24 CM	52224-79A	52224-79B
	33 CM	52255-69A	52255-69B
	36 CM	52236-69A	52236-69B
	45 CM	52155-69A	52155-69B
	24 CM	52224-69A	52224-69B
	33 CM	52255-70A	52255-70B
	36 CM	52236-70A	52236-70B
	45 CM	52155-70A	52155-70B
	24 CM	52224-70A	52224-70B
	33 CM	52255-68A	52255-68B
	36 CM	52236-68A	52236-68B
	45 CM	52155-68A	52155-68B
	24 CM	52224-68A	52224-68B
	33 CM	52255-71A	52255-71B
	36 CM	52236-71A	52236-71B
	45 CM	52155-71A	52155-71B
	24 CM	52224-71A	52224-71B
	33 CM	52255-72A	52255-72B
	36 CM	52236-72A	52236-72B
	45 CM	52155-72A	52155-72B
	24 CM	52224-72A	52224-72B
	33 CM	52255-73A	52255-73B
	36 CM	52236-73A	52236-73B
	45 CM	52155-73A	52155-73B
	24 CM	52224-73A	52224-73B



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

10 MM Ø Dissectors & Graspers

Rigid One-Piece with Flush Port Rotatable One-Piece with Flush Port

	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
	33 CM	52102-02A	52102-02B
	36 CM	52103-02A	52103-02B
	45 CM	52101-02A	52101-02B
10MM, MARYLAND GRASPER/DISSECTOR	33 CM	52102-06A	52102-06B
	36 CM	52103-06A	52103-06B
	45 CM	52101-06A	52101-06B
10MM, MIXTER GRASPER/DISSECTOR, RIGHT ANGLE	33 CM	52102-07A	52102-07B
	36 CM	52103-07A	52103-07B
	45 CM	52101-07A	52101-07B
10MM, MIXTER GRASPER/DISSECTOR, SHORT RIGHT ANGLE	33 CM	52102-59A	52102-59B
	36 CM	52103-59A	52103-59B
	45 CM	52101-59A	52101-59B
10MM, MIXTER GRASPER/DISSECTOR, 45° ANGLE	33 CM	52102-08A	52102-08B
	36 CM	52103-08A	52103-08B
	45 CM	52101-08A	52101-08B
10MM, BABCOCK GRASPER, FENESTRATED TIP	33 CM	52102-16A	52102-16B
	36 CM	52103-16A	52103-16B
	45 CM	52101-16A	52101-16B
US BABCOCK GRASPER, SPRING TENSION, ATRAUMATIC	33 CM	52102-17A	52102-17B
	36 CM	52103-17A	52103-17B
	45 CM	52101-17A	52101-17B
10MM, ALLIS GRASPER, FINE AXIAL SERRATIONS	33 CM	52102-21A	52102-21B
	36 CM	52103-21A	52103-21B
	45 CM	52101-21A	52101-21B
10MM, DEBAKEY GRASPER, 40 MM JAW	33 CM	52102-22A	52102-22B
	36 CM	52103-22A	52103-22B
	45 CM	52101-22A	52101-22B
10MM, DEBAKEY GRASPER, 60 MM JAW	33 CM	52102-23A	52102-23B
	36 CM	52103-23A	52103-23B
	45 CM	52101-23A	52101-23B



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

10 MM Ø Dissectors & Graspers

		Rigid One-Piece with Flush Port	Rotatable One-Piece with Flush Port
	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
	33 CM	52102-24A	52102-24B
	36 CM	52103-24A	52103-24B
	45 CM	52101-24A	52101-24B
10MM, KOCHER GRASPER, 60MM JAW, AXIAL SERRATIONS	33 CM	52102-25A	52102-25B
	36 CM	52103-25A	52103-25B
	45 CM	52101-25A	52101-25B
10MM, DUVAL GRASPER, 40 MM JAW, AXIAL SERRATIONS	33 CM	52102-26A	52102-26B
	36 CM	52103-26A	52103-26B
	45 CM	52101-26A	52101-26B
10MM, CLAW FORCEPS, 2 X 3 TEETH	33 CM	52102-36A	52102-36B
	36 CM	52103-36A	52103-36B
	45 CM	52101-36A	52101-36B
10MM, SPOON FORCEPS, STANDARD TIP	33 CM	52102-37A	52102-37B
	36 CM	52103-37A	52103-37B
	45 CM	52101-37A	52101-37B
10MM, SPOON FORCEPS, LONG TIP	33 CM	52102-38A	52102-38B
	36 CM	52103-38A	52103-38B
	45 CM	52101-38A	52101-38B
10MM, SCHROEDER TENACULUM FORCEPS, SINGLE TOOTH	33 CM	52102-45A	52102-45B
	36 CM	52103-45A	52103-45B
	45 CM	52101-45A	52101-45B
10MM, SCHROEDER TENACULUM FORCEPS "WISAP STYLE," SINGLE TOOTH W/ EXTRA GUIDE TOOTH	33 CM	52102-94A	52102-94B
	36 CM	52103-94A	52103-94B
	45 CM	52101-94A	52101-94B
10MM, STRONG ATRAUMATIC CLAMP, 40MM STRAIGHT JAW	33 CM	52102-46A	52102-46B
	36 CM	52103-46A	52103-46B
	45 CM	52101-46A	52101-46B
10MM, STRONG ATRAUMATIC CLAMP, 40MM CURVED JAW	33 CM	52102-47A	52102-47B
	36 CM	52136-47A	52136-47B
	45 CM	52101-47A	52101-47B



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

10 MM Forceps, Biopsys, Tenaculums, and Scissors

		Rigid One-Piece with Flush Port	Rotatable One-Piece with Flush Port
	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
	33 CM	52102-24A	52102-24B
	36 CM	52103-24A	52103-24B
	45 CM	52101-24A	52101-24B
	33 CM	52102-25A	52102-25B
	36 CM	52103-25A	52103-25B
	45 CM	52101-25A	52101-25B
	33 CM	52102-26A	52102-26B
	36 CM	52103-26A	52103-26B
	45 CM	52101-26A	52101-26B
	33 CM	52102-36A	52102-36B
	36 CM	52103-36A	52103-36B
	45 CM	52101-36A	52101-36B
	33 CM	52102-37A	52102-37B
	36 CM	52103-37A	52103-37B
	45 CM	52101-37A	52101-37B
	33 CM	52102-38A	52102-38B
	36 CM	52103-38A	52103-38B
	45 CM	52101-38A	52101-38B
	33 CM	52102-45A	52102-45B
	36 CM	52103-45A	52103-45B
	45 CM	52101-45A	52101-45B
	33 CM	52102-94A	52102-94B
	36 CM	52103-94A	52103-94B
	45 CM	52101-94A	52101-94B
	33 CM	52102-46A	52102-46B
	36 CM	52103-46A	52103-46B
	45 CM	52101-46A	52101-46B
	33 CM	52102-47A	52102-47B
	36 CM	52136-47A	52136-47B
	45 CM	52101-47A	52101-47B
	33 CM	52102-64A	52102-64B
	36 CM	52136-64A	52136-64B
	45 CM	52101-64A	52101-64B
	33 CM	52102-67A	52102-67B
	36 CM	52136-67A	52136-67B
	45 CM	52101-67A	52101-67B



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Additional Handle Styles



TYPE 1:
INSULATED WITH RF-
POST



TYPE 2:
INSULATED W/ RF-
POST & RATCHET



TYPE 3:
NON-INSULATED
STAINLESS STEEL



TYPE 4:
NON-INSULATED
STAINLESS STEEL W/
RATCHET



TYPE 5:
LIGHTWEIGHT
CARBON-FIBER W/ RF
POST



TYPE 6:
LIGHTWEIGHT
CARBON-FIBER
w/ RF-POST &
RATCHET



TYPE 7:
LIGHTWEIGHT
CARBON-FIBER
W/ RF POST,
RATCHET /
RATCHET
DEFEATING



TYPE 8:
LIGHTWEIGHT CARBON
-FIBER W/ RF POST,
FRENCH ON/OFF
RATCHET



TYPE 9:
LIGHTWEIGHT
CARBON-FIBER W/ RF
POST, W/ LARGE
FINGER LOOPS



TYPE 10:
LIGHTWEIGHT
CARBON-FIBER W/ RF
POST, RATCHET /
RATCHET DEFEATING
W/ LARGE FINGER
LOOPS



TYPE 11:
STAINLESS STEEL RING
AXIAL HANDLE



TYPE 12:
STAINLESS STEEL RING
HANDLE W/ RATCHET



TYPE 13:
STAINLESS STEEL AXIAL HANDLE
W/ PRECISE ENGAGEMENT
RATCHET



TYPE 14:
INSULATED STAINLESS
STEEL LINVATEC-WECK
STYLE W/ LONG
SHANKS



TYPE 15:
NON-INSULATED
STAINLESS STEEL
LINVATEC-WECK STYLE
W/ QUICK RELEASE
RATCHET & LONG
SHANKS



TYPE 16:
NON-INSULATED
STAINLESS STEEL W/
SPRING ASSIST
CLOSING



Type 17:
Lightweight Carbon
-Fiber w/ RF Post



TYPE 18:
CARBON-FIBER W/ RF
POST, RATCHET /
RATCHET DEFEATING