



MILLENNIUM
SURGICAL

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Bipolar Bayonet Forceps Non-Stick Tips

Our insulated forceps are made in the U.S.

Compare to Stryker's Silverglide*, Codman's Mallis*, Integra's Buzz Lite*!

Malis-Style Bayonet Bipolar Forceps



Item Number	Description	Length	Tip Size
100-S03648NS	Malis-style Bayonet	7.5" (19 cm)	.25 mm
100-S03651NS	Malis-style Bayonet	7.5" (19 cm)	.5 mm
100-S03653NS	Malis-style Bayonet	7.75" (19.7 cm)	1.0 mm
100-S03655NS	Malis-style Bayonet	7.5" (19 cm)	1.5 mm

S02

S03



Item Number	Description	Length	Tip Size
100-S03674TR	Malis-style Bayonet	8.5" (21.8 cm)	.25 mm
100-S03675TR	Malis-style Bayonet	8.5" (21.8 cm)	.5 mm
100-S03676TF	Malis-style Bayonet	8.5" (21.8 cm)	.7 mm
100-S03677TR	Malis-style Bayonet	8.5" (21.8 cm)	1.0 mm
100-S03678TR	Malis-style Bayonet	8.5" (21.8 cm)	1.5 mm

S02


S03



Item Number	Description	Length	Tip Size
100-S03680NS	Malis-style Bayonet	8.75" (22.4 cm)	.25 mm
100-S03681NS	Malis-style Bayonet	8.75" (22.4 cm)	.5 mm
100-S03682NS	Malis-style Bayonet	8.75" (22.4 cm)	.7 mm
100-S03683NS	Malis-style Bayonet	8.75" (22.4 cm)	1.0 mm
100-S03685NS	Malis-style Bayonet	8.75" (22.4 cm)	1.5 mm

S02

S03



Item Number	Description	Length	Tip Size
100-S03695NS	Malis-style Bayonet	10.75" (27.3 cm)	.25 mm
100-S03696NS	Malis-style Bayonet	10.75" (27.3 cm)	.5 mm
100-S03697NS	Malis-style Bayonet	10.75" (27.3 cm)	.7 mm
100-S03698NS	Malis-style Bayonet	10.75" (27.3 cm)	1.0 mm
100-S03699NS	Malis-style Bayonet	10.75" (27.3 cm)	1.5 mm

S02


S03



Symbols



Tip Sizing Chart





MILLENNIUM
SURGICAL

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Bipolar Bayonet Forceps Non-Stick Tips

Key and Keyfin Bayonet Bipolar Forceps

Available ADD-ON



S03

Item Number	Description	Length	Tip Size
100-S03230FG	Keyfin Bayonet	8.00" (20.3 cm)	.4 mm
100-S03233FG	Keyfin Bayonet	8.00" (20.3 cm)	.7 mm
100-S03235FG	Keyfin Bayonet	8.00" (20.3 cm)	1.0 mm
100-S03237FG	Keyfin Bayonet	8.00" (20.3 cm)	1.2 mm
100-S03238FG	Keyfin Bayonet	8.00" (20.3 cm)	1.5 mm

NS FG



S03

Item Number	Description	Length	Tip Size
100-S03240FG	Keyfin Bayonet	9.00" (22.7 cm)	.4 mm
100-S03243FG	Keyfin Bayonet	9.00" (22.7 cm)	.7 mm
100-S03245FG	Keyfin Bayonet	9.00" (22.7 cm)	1.0 mm
100-S03247FG	Keyfin Bayonet	9.00" (22.7 cm)	1.2 mm
100-S03248FG	Keyfin Bayonet	9.00" (22.7 cm)	1.5 mm

NS FG



S03

Item Number	Description	Length	Tip Size
100-S03260FG	Keyfin Bayonet	10.00" (25.4 cm)	.4 mm
100-S03263FG	Keyfin Bayonet	10.00" (25.4 cm)	.7 mm
100-S03265FG	Keyfin Bayonet	10.00" (25.4 cm)	1.0 mm
100-S03267FG	Keyfin Bayonet	10.00" (25.4 cm)	1.2 mm
100-S03268FG	Keyfin Bayonet	10.00" (25.4 cm)	1.5 mm

NS FG



S03

Item Number	Description	Length	Tip Size
100-S03280FG	MIS Keyfin Bayonet	11.00" (27.9 cm)	.4 mm
100-S03281FG	MIS Keyfin Bayonet	11.00" (27.9 cm)	.7 mm
100-S03285FG	MIS Keyfin Bayonet	11.00" (27.9 cm)	1.0 mm
100-S03286FG	MIS Keyfin Bayonet	11.00" (27.9 cm)	1.2 mm
100-S03288FG	MIS Keyfin Bayonet	11.00" (27.9 cm)	1.5 mm

NS FG

Symbols

S02 S03
⊗ NS SS S T R SP FG
 Non-insulated Insulated Stainless tip Non-stick tip 2 Stops 1 Stop Titanium Irrigation Stop Pin Fin Guide

Tip Sizing Chart

