



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Reusable Cables

For monopolar and bipolar instruments and generators

- Repairable and Reusable - autoclavable and for use with laparoscopic and electrosurgical instruments.
- Jacketing of the cable is made of high-grade silicone.
- Machined (not molded) connecting plugs on both ends of the cables for strength and repairability.
- Manufactured to original manufacturer's standards and compatible with electrosurgical units.



GENERATOR END



GENERATOR END

INSTRUMENT END

LIST NO. 5021	4MM DIA. PLUG (UNIT) & 4MM DIA. INST. PLUG		LIST NO. 5009	8MM DIA. "BOVIE" PLUG (UNIT) & 4MM DIA. INST. PLUG
LIST NO. 5022	4MM DIA. PLUG (UNIT) & OLSEN TYPE SCREW- ON INST. PLUG		LIST NO. 5023	8MM DIA. "BOVIE" PLUG (UNIT) & OLSEN TYPE SCREW-ON INST. PLUG
LIST NO. 5029	4MM DIA. PLUG (UNIT) & OLYMPUS TYPE INST. PLUG		LIST NO. 5030	8MM DIA. "BOVIE" PLUG (UNIT) & OLYMPUS TYPE INST. PLUG
LIST NO. 5260	4MM DIA. PLUG (UNIT) & ACMI TYPE INST. PLUG		LIST NO. 5261	8MM DIA. "BOVIE" PLUG (UNIT) & ACMI TYPE INST. PLUG
LIST NO. 5268	4MM DIA. PLUG (UNIT) & TAPERED "DOWNS" TYPE INST. PLUG		LIST NO. 5269	8MM DIA. "BOVIE" PLUG (UNIT) & TAPERED "DOWNS" TYPE PLUG
LIST NO. 52005-29	4MM DIA. PLUG (UNIT) & 3MM DIA. FEMALE INST. PLUG		LIST NO. 52005-30	8MM DIA. "BOVIE" PLUG (UNIT) & 3MM DIA. FEMALE INST. PLUG
LIST NO. 52005-31	4MM DIA. PLUG (UNIT) & 4MM DIA. FEMALE INST. PLUG		LIST NO. 52005-32	8MM DIA. "BOVIE" PLUG (UNIT) & 4MM DIA. FEMALE INST. PLUG
LIST NO. 52005-35	4MM DIA. PLUG (UNIT) & CORSON NEEDLE ELECTRODE INST. PLU		LIST NO. 52005-36	8MM DIA. "BOVIE" PLUG (UNIT) & CORSON NEEDLE ELECTRODE INST. PLUG
LIST NO. 52005-37	4MM DIA. PLUG (UNIT) & DUAL STAGE MALE (STORZ) INST. PLUG		LIST NO. 52005-38	8MM DIA. "BOVIE" PLUG (UNIT) & DUAL STAGE MALE (STORZ) INST. PLUG
LIST NO. 52005-47	4MM DIA. PLUG (UNIT) & RESECTOSCOPE INST. PLUG		LIST NO. 52005-48	8MM DIA. "BOVIE" PLUG (UNIT) & RESECTOSCOPE INST. PLUG
LIST NO. 52005-49	4MM DIA. PLUG (UNIT) & 4MM DIA. RECESSED PIN (WOLF/CORSON/STORZ) INST. PLUG		LIST NO. 52005-50	8MM DIA. "BOVIE" PLUG (UNIT) & 4MM DIA. RECESSED PIN (WOLF/CORSON/STORZ) INST.
LIST NO. 52005-69	4MM DIA. PLUG (UNIT) & 4MM DIA. SELF- RETAINING PLUG FOR STORZ SOLID POST INST.		LIST NO. 52005-70	8MM DIA. "BOVIE" PLUG (UNIT) & 4MM DIA. SELF- RETAINING PLUG FOR STORZ SOLID POST INST.
LIST NO. 52005-71	4MM DIA. PLUG (UNIT) & 4MM DIA. INSULATED OUTER SLEEVE PLUG FOR V.MUELLER INST		LIST NO. 52005-72	8MM "BOVIE" PLUG (UNIT) & 4MM INSULATED OUTER SLEEVE PLUG FOR V.MUELLER INST.
LIST NO. 52005-73	4MM DIA. EXTRA-LONG PLUG (UNIT) & 4MM DIA. EXTRA-LONG PLUG: STORZ INST.		LIST NO. 52005-74	8MM DIA. "BOVIE" PLUG (UNIT) & 4MM DIA. EXTRA-LONG PLUG: STORZ INST
LIST NO. 52005-75	4MM DIA. PLUG (UNIT) & 4MM DIA. RIGHT ANGLE FEMALE INST. PLU .		LIST NO. 52005-76	8MM DIA. "BOVIE" PLUG (UNIT) & 4MM DIA. RIGHT ANGLE FEMALE INST. PLUG
LIST NO. 52005-77	4MM DIA. PLUG (UNIT) & MALE LUER LOCK INST. PLUG		LIST NO. 52005-78	8MM DIA. "BOVIE" PLUG (UNIT) & MALE LUER LOCK INST. PLUG



LIST NO.
5414

ADAPTER FOR MALE "BOVIE"
GENERATOR PLUG & 4MM FE-
MALE BANANA CABLE PLUG



LIST NO.
5415

ADAPTER FOR MALE "BOVIE"
GENERATOR PLUG & DISPOS-
ABLE INSTRUMENT PLUG



LIST NO.
5416

ADAPTER FOR MALE
BANANA GENERATOR PLUG &
DISPOSABLE INSTRUMENT PLUG



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Reusable Cables

For monopolar and bipolar instruments and generators



GENERATOR END

INSTRUMENT END



GENERATOR END

LIST NO. 52005-21 / 5302	ELMED/WOLF PLUG (UNIT) & ELMED (DIN) INSTRUMENT PLUG		LIST NO. 52005-22 / 5230	DUAL-BANANA PLUG (UNIT) & ELMED (DIN) INSTRUMENT PLUG
LIST NO. 52005-19	ELMED/WOLF PLUG (UNIT) & ELMED TYPE III INST. PLUG		LIST NO. 52005-20	DUAL-BANANA PLUG (UNIT) & ELMED TYPE III INST. PLUG
LIST NO. 52005-23	ELMED/WOLF UNIT PLUG & WOLF INST. PLUG		LIST NO. 52005-24	DUAL-BANANA PLUG (UNIT) & WOLF INST. PLUG
LIST NO. 52005-25	ELMED/WOLF PLUG (UNIT) & STORZ/WISAP/JARIT METAL INST. PLUG		LIST NO. 52005-26	DUAL-BANANA PLUG (UNIT) & STORZ/WISAP/JARIT METAL PLUG
LIST NO. 52005-17	ELMED/WOLF PLUG (UNIT) & STORZ/WISAP BARREL STYLE INST. PLUG		LIST NO. 52005-18	DUAL-BANANA PLUG (UNIT) & STORZ/WISAP BARREL STYLE INST. PLUG
LIST NO. 52005-39	ELMED/WOLF PLUG (UNIT) & KLI/CIRCON INST. PLUG		LIST NO. 52005-40	DUAL-BANANA PLUG (UNIT) & KLI/CIRCON INST. PLUG
LIST NO. 52005-41	ELMED/WOLF PLUG (UNIT) & EVEREST/COBALT INST. PLUG		LIST NO. 52005-42	DUAL-BANANA PLUG (UNIT) & EVEREST/COBALT INST. PLUG
LIST NO. 52005-43	ELMED/WOLF PLUG (UNIT) & DAVOL INST. PLUG		LIST NO. 52005-44	DUAL-BANANA PLUG (UNIT) & DAVOL INST. PLUG
LIST NO. 52005-45 / 5991	ELMED/WOLF PLUG (UNIT) & AESCULAP/DAVOL/MARINA 2-PIN INST. PLUG		LIST NO. 52005-46 / 5990	DUAL-BANANA PLUG (UNIT) & AESCULAP/DAVOL/MARINA 2-PIN INST. PLUG
LIST NO. 52005-60	ELMED/WOLF PLUG (UNIT) & LINVATEC/WECK/PILLING INST. PLUG		LIST NO. 52005-61	DUAL-BANANA PLUG (UNIT) & LINVATEC/WECK/PILLING INST. PLUG
LIST NO. 52005-91	ELMED/WOLF PLUG (UNIT) & JARIT, ERBE INST. PLUG		LIST NO. 52005-92	DUAL-BANANA PLUG (UNIT) & JARIT, ERBE INST. PLUG
LIST NO. 52005-13	ELMED/WOLF PLUG (UNIT) & ELMED BI-CUT SCISSORS		LIST NO. 52005-14	DUAL-BANANA PLUG (UNIT) & ELMED BI-CUT SCISSORS

BIPOLAR CABLE ADAPTERS



LIST NO.
5352

ADAPTER FOR ELMED / WOLF
GENERATOR PLUG & 4MM
FEMALE BANANA CABLE PLUGS



LIST NO.
52005-83

ADAPTER FOR DUAL MALE BANANA
GENERATOR PLUGS & ELMED /
WOLF RECEPTACLE PLUG



LIST NO.
52005-84

ADAPTER FOR DUAL MALE
BANANA GENERATOR PLUGS &
WOLF KLEPPINGER FORCEPS PLUG