



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

I/A Handpieces and Tips

92-IA21



Reduce The Risk of TASS

The handle disassembles for thorough cleaning and inspection and is designed with fewer nooks and crannies, reducing the risk of contamination. Manufactured of corrosion-free titanium.

Improve Procedural Outcome

Avoid glare! Handle and tips are made of light-weight, glare-free titanium. A large selection of tip styles and port sizes are available for use in standard and micro cataract surgeries. Tips can be rotated intraoperatively!

- Fully compatible with all phaco machines and most brands of existing I/A Handpiece systems.
- Compatible with these phaco systems: Alcon*, B&L Storz*, AMO*, Abbot*, Starr*
- Compatible with these handpiece systems: B&L Storz*, AMO*, MST*, Ambler*
- Made in America!

Part Number	Product	Description
92-IA21	I/A Handpiece Handle	Disassembles for easy cleaning.
92-IA22	Straight tip	.3mm port. For initial cortical removal.
92-IA225	Straight tip	.5mm port. For initial cortical removal.
92-IA24	45 Degree Tip	.3mm port. For removing periperal cortical material.
92-IA25	90 Degree Tip	.3mm port. For removing periperal cortical material.
92-IA26	150 Degree Tip	.3mm port. For aspiration of 12 o'clock cortex.
92-IA23	Curved Tip	.3 mm port. For deep-set eyes.
92-IA235	Curved Tip	.5 mm port. For deep-set eyes.
92-IA27	Straight Tip for Silicone Sleeve	.300. For small incision cortical cleanup.
92-IA29	45 Degree Tip for Silicone Sleeve	Useful for peripheral cortical removal.
92-IA28	CVD Tip for Silicone Sleeve	.300. Useful for deep-set eyes.