



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

I/A Handpieces and Tips

92-IA21



Improve Procedural Outcome

Avoid glare! Handle and tips are made of light-weight, glare-free titanium. A large selection of tip styles and port sizes are available for use in standard and micro cataract surgeries. Tips can be rotated intraoperatively!



- Fully compatible with all phaco machines and most brands of existing I/A Handpiece systems.
- Compatible with these phaco systems: Alcon*, B&L Storz*, AMO*, Abbot*, Starr*
- Compatible with these handpiece systems: B&L Storz*, AMO*, MST*, Ambler*
- Made in America!

Part Number	Product	Description
92-IA21	I/A Handpiece Handle	Disassembles for easy cleaning.
92-IA22	Straight tip	.3mm port. For initial cortical removal.
92-IA225	Straight tip	.5mm port. For initial cortical removal.
92-IA24	45 Degree Tip	.3mm port. For removing periperal cortical material.
92-IA25	90 Degree Tip	.3mm port. For removing periperal cortical material.
92-IA26	150 Degree Tip	.3mm port. For aspiration of 12 o'clock cortex.
92-IA23	Curved Tip	.3 mm port. For deep-set eyes.
92-IA235	Curved Tip	.5 mm port. For deep-set eyes.
92-IA27	Straight Tip for Silicone Sleeve	.300. For small incision cortical cleanup.
92-IA29	45 Degree Tip for Silicone Sleeve	Useful for peripheral cortical removal.
92-IA28	CVD Tip for Silicone Sleeve	.300. Useful for deep-set eyes.

*The above referenced company names are registered trademarks of the respective owners of the companies. Millennium Surgical Corp. and surgicalinstruments.com are not affiliated with the above companies. Millennium Surgical does not buy or distribute products from the above companies.