



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Backflush Instruments

Reusable Backflush Instruments



5146 - Titanium Backflush Instrument - Passive
1 per pack, non-sterile



5153 - Titanium Backflush Instrument - Active
1 per pack, non-sterile

Disposable Backflush Instruments



5197 - Disposable Backflush Instrument
Passive. 5 per box, sterile



5196 - Disposable Backflush Instrument
Active. 5 per box, sterile

Reservoirs for Backflush Instruments - Sterile



5147 - Passive Replacement Reservoir for B/F Instrument - 10 / box



5159 - Active Replacement Reservoir for B/F Instrument - 10 / box

Silicone Tip Cannulas

20Ga



5017 - Brush Tip Cannula, 20Ga - 0.9mm silicone brush tip


5063 - Soft Tip Cannula, 20Ga - 2.0mm silicone tip

5193 - Diamond Dusted Soft Tip Cannula, 20Ga - 1.0mm silicone tip
(for triamcinolone removal)

5151 - Charles Flute Cannula, 20Ga - Blunt tip



Brush Tip Cannula



Soft Tip Cannula



Diamond Dusted Cannula



Charles Flute Cannula

23Ga

5162 - Brush Tip Cannula, 23Ga - 0.9mm silicone brush tip

5161 - Soft Tip Cannula, 23Ga - 2.0mm silicone tip

5192 - Diamond Dusted Soft Tip Cannula, 23Ga - 1.0mm silicone tip
(for triamcinolone removal)

5164 - Charles Flute Cannula, 23Ga - Blunt tip

25Ga

5160 - Brush Tip Cannula, 25Ga - 0.6mm silicone brush tip

5152 - Soft Tip Cannula, 25Ga - 2.0mm silicone tip

5156 - Charles Flute Cannula, 25Ga - Blunt tip



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Retinal Cannulas

Infusion Cannulas, 20Ga

- 5066 - Infusion Cannula 20Ga, 2.5mm - 5 per box
- 5010 - Infusion Cannula 20Ga, 4.0mm - 5 per box
- 5011 - Infusion Cannula 20Ga, 6.0mm - 5 per box




Silicone Oil Cannulas, 20Ga

- 5012-Silicone Oil Infusion Cannula 20Ga, 4.0mm-5 / box
- 5013-Silicone Oil Infusion Cannula 20Ga, 6.0mm-5 / box



Self-Retaining Sutureless Cannulas 20Ga /23Ga

- Sutureless fixation
- Retaining ridge holds cannula in place
- Stainless steel cannula hub



- 5116-Self-Retaining Infusion Cannula 20Ga, 4.0mm - 5 per box
- 5144-Self-Retaining Infusion Cannula 23Ga, 4.0mm - 5 per box
- 5219-Self-Retaining Infusion Cannula 20Ga, 6.0mm - 5 per box
- 5221-Self-Retaining Infusion Cannula 23Ga, 6.0mm - 5 per box


Self-Retaining Silicone Oil Sutureless Cannulas

- 5117-Self-Retaining Silicone Oil Cannula 20Ga, 4.0mm - 5/box
- 5165-Self-Retaining Silicone Oil Cannula 23Ga, 4.0mm - 5/box
- 5220-Self-Retaining Silicone Oil Cannula 20Ga, 6.0mm - 5/box
- 5222-Self-Retaining Silicone Oil Cannula 23Ga, 6.0mm - 5/box



Dual Bore PFC Cannulas 20Ga and 23Ga

- 5191 - Dual Bore PFC Cannula, 20Ga - 5 per box
- 5203 - Dual Bore PFC Cannula, 23Ga - 5 per box
- Simultaneous infusion of heavy liquids and aspiration of intraocular fluids
- Dual bore enables to control constant intraocular pressure during injection



Dual Bore Subretinal BSS Injection Needle, 20Ga

- 5194-Dual Bore BSS Injection Needle 0.1mm / 41Ga tip
- 1 per box
- Enables to control subretinal injection of BSS
- Dual bore cannula combined with aspiration capability

