




**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM


UNMATCHED SERVICE · UNLIMITED INSTRUMENTS








822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Diamond-Like Sapphire Knives for Eye Surgery



- Crisp, sharp & precise incision for faster healing, less scar
- Enhanced performance due to high density of gemstone
- Less damage, no metal particles left in tissue
- Lightweight & durable Titanium handle



Code	Description	B	A	Specifications				
				D	C	L	E	F
	60-FSK01A Sideport	1.0	15°	0.20	24			
	60-FSK01B Sideport	1.0	30°	0.20	24			
	60-FSK02 Two sided Keratome	2.80 mm	60°	0.30	2.42	4.50		
	60-FSK03 Two sided Keratome	3.00 mm	60°	0.30	2.60	4.50		
	60-FSK04 Two sided Keratome	3.20 mm	60°	0.30	2.77	4.50		
	60-FSK05A Crescent	2.00 mm		0.25				
	60-FSK05B Crescent	2.80 mm		0.25				
	60-FSK06A Trapezoid	2.90 mm	2.10	0.25	3.00	7.50	5.95	
	60-FSK06B Trapezoid	3.20 mm	2.35	0.25	3.00	7.50	5.95	0.48
	60-FSK06C Trapezoid	3.40 mm	2.50	0.25	3.00	7.50	5.95	0.48
	60-FSK07 MVR							
	60-FSK08 Trifacet	1.00	30°	0.20		7.50	45°	
	60-FSK09 All sided Keratome	2.80 mm	60°			7.00	2.77	
	60-FSK10 All sided Keratome	3.00 mm	60°			7.00	2.77	
	60-FSK11 All sided Keratome	3.20 mm	60°			7.00	2.77	