



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Lap Needle Electrodes

Monopolar Needle Electrodes for Laparoscopic Surgery
Exclusive Retractable Needle Electrodes for Precise Laparoscopic
Cutting and Coagulation



FIGURE 1



FIGURE 2

Item Number	Description
72-5220553	30 CM Insulated cutting electrode, retractable needle, according to Dr. Valle, 5 MM diameter shaft with 3 MM tip for better vision, set screw adjustment (figure 1).
72-5220554	Replacement Stylet for 72-5220553
72-5220569	30 CM Insulated cutting electrode, retractable needle, with irrigation by trumpet valve, according to Dr. Valle, 5 MM diameter shaft, set screw adjustment.
72-5220570	Irrigation sleeve only, for 72-5220569
72-52205110	30 CM Insulated cutting electrode, retractable needle, 5 MM diameter shaft with 3 MM tip for better vision, spring loaded push button (figure 2).



FIGURE 3

Item Number	Description
72-5220548	30 CM Insulated cutting electrode, with chuck for micro needles, according to Dr. Valle, 5 MM diameter shaft (figure 3).
72-5200533	Replacement 150 Micron Needle for 72-5220548
72-5200534	Replacement 500 Micron Needle for 72-5220548

The same monopolar needle electrodes are available in 44 CM.
Please request order number and price.