



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Xenon Lightsources



System Warranty

Light Source = 3 Years Cable = 12 Months
Headlight = 18 Months Bulb = 700 Hours



White Star 175 Watt Xenon Light Source (Above)

Also Available:
Dual White Star 175 Watt Xenon Light Source
(Right, Shown with Universal Floorstand)



White Sun® 300 Watt Xenon Light Source (Above)

Also Available:
Dual White Sun® 300 Watt Xenon Light Source
(Right, Shown with Universal Floorstand)



White Sun® 300 Watt Xenon Light Source on white stand



PFM Light Sources feature fast, easy, on-the-spot lamp replacement

Call us today to request a demo unit!
1-800-600-0428

Order No.	Description
XL 300 N	White Sun® 300 Watt Xenon Light Source
XL 300 N DUAL	Dual White Sun® 300 Watt Xenon Light Source
XL 175 N	White Star 175 Watt Xenon Light Source
XL 175 N DUAL	Dual White Star 175 Watt Xenon Light Source
XL 200 (not shown)	Small Sun 200 Watt Quartz Halogen Light Source
X 3114	Blue Universal Floor Stand
X 3114 W	White / Chrome Universal Floor Stand



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Xenon Head Lights



White Sun Series

White Sun® Headlight

White Sun® HL w/ Sweatband



White Sun Series



Minnie Series

Minnie Headlight

Minnie HL w/ Sweatband



Minnie Series

Ordering Guide

White Sun Series		Traditional Series		Minnie Series	
X 781 HL	Headlight with Rigid Headband	Single Cable (not shown)		X 779 HL	Headlight with Rigid Headband
X 781 HL MOD	Headlight with Sweatband	X 772 HL	Headlight with Rigid Headband	X 779 HL MOD	Headlight with Sweatband
PEM 3095 BIF	3mm × 9.5' Bifurcated Cable	X 770 5	5mm × 7.5' Single Cable	X 778 1	Sweatbands, 25/bx.
PEM 3075 BIF	3mm × 7.5' Bifurcated Cable	X 775 1	5mm × 10' Single Cable	<div style="border: 1px solid black; padding: 2px; display: inline-block;"> Each Minnie Series headlight comes with a 9' cable attached. </div>	
X 778 3	Sweatbands, 25/bx.	X 3075	3mm × 7.5' Single Cable		
White Sun Documentation		X 3010	3mm × 10' Single Cable		
X 775 HL	White Sun Headlight with video (not shown)	Bifurcated Cable (right)			
Portable Series		X 773 HL	Headlight with Rigid Headband		
X 777 HL	Portable battery operated Headlight (not shown)	X 773 2	5mm × 7.5' Bifurcated Cable		
		X 773 3	5mm × 10' Bifurcated Cable		
		3075 BIF	3mm × 7.5' Bifurcated Cable		
		3095 BIF	3mm × 9.5' Bifurcated Cable		



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Illuminarc® Portable Headlamp



- Freedom of movement around the operating room
- Bright light intensity for all surgical procedures
- Cost savings compared to plug-in headlight systems
- Illuminarc® light technology gives true 5500k light in the body cavity as compared to the blue light of most other portable LED light technology

The XL 250 CS includes:

- 1 - Headlight
- 1 - Fiber Optic Cable
- 1 - Light Source
- 1 - Belt
- 1 - Backpack
- 1 - Battery Charger
- 2 - 2 Hour Batteries
- Or
- 2 - 3 Hour batteries





**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Xenon Head Lights

XL 150 CS LED Headlight System



This LED headlight produces an unparalleled bright, white light with superior battery life so you can work without interruption.

The Cool Air design Head Frame will create a new wearing experience - lightweight, comfortable, flexible and durable.

Luminescent and feather weight 10W LED lighting system. When used at full light intensity the battery will last 2 hours 40 minutes, providing you with 34,000 lux of light output to the surgical site.

LED Headlight Charger A compact and highly efficient charger will fully charge an empty fuel cell battery in 3 ½ hours and stop charging once fully charged. Accepts 100-240V AC, 50/60 Hz input for worldwide use.

10w LED Producing Bright, White Light

Color temperature: 4500°K

Typical LED life: >50,000 hours

Newly Designed Lightweight Headband with Rear Cranial Support for added comfort

Headlight Module with 70mm Fixed Spot at 16"

Variable Intensity Control: 0 – 34,000 Lux at 16"

Headband with LED Module Weight – 9 oz

Battery and Holster Weight – 12 oz

Battery Life – 2 hrs 40 min @ Full Intensity, 5 hrs 30 min at Half Intensity

Battery Charger Accepts 100-240V AC, 50/60 Hz input for worldwide operation

Battery Recharge Time – 3.5 hrs

Warranty Three (3)Years on LED Headlight and Six (6) Months on Battery

System Includes

- 1 LED Headlight With Head Frame
- 1 Clip-on Power Pack with rechargeable Battery
- 1 2" Heavyweight Polypropylene Quick Release Belt
- 1 Single Bay Battery Charger
- 2 Garment Clips



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE • UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Xenon Head Lights

Item Numbers:

X 781 HL NWR	White Ray 300 Watt Xenon Light Source
PFM 3095 BIF	10 ft. Bifurcated Fiber Optic Cable
XL 300 NWR	White Ray Headlight
X 3114W NWR	Floor Stand
X 7787 NWR	Storage Bag
XL 300 CS NWR	Complete System

Warranty:

Light Source:	4 Years
Xenon Lamp:	1000 Hours
Headlight:	18 Months
Fiber Optic Cable:	12 Months

