



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Muscle Biopsy Needles


UCH Type Skeletal Muscle Biopsy Needles - Stainless Steel

A: Outer cannula with printed tip and fenestrated shaft with dual ring handle (order by shaft Ø).

B: Inner cannula with cutting tip and female luer fitting for suction application (finger control suction) and ring handle.

C: Blunt trocar for removal of sample

D: Sealing cap for inner cannula (optional use)



Product Number	Description	Standard Wire Gauge
72-238064	UCH Needle 3mm Ø x 120mm	10 swg x 4.75"
72-238065	UCH Needle 4mm Ø x 120mm	8 swg x 4.75"
72-238066	UCH Needle 4mm Ø x 80mm	8 swg x 3.25"
72-238067	UCH Needle 4.5mm Ø x 120mm	7 swg x 4.75"
72-238068	UCH Needle 5mm Ø x 120mm	6 swg x 4.75"
72-238069	UCH Needle 6mm Ø x 120mm	5 swg x 4.75"


Bergstrom Style Skeletal Muscle Biopsy Needles - Stainless Steel

A: Outer cannula with printed tip and fenestrated shaft with dual ring handle (order by shaft Ø).

B: Inner cannula with cutting tip and female luer fitting for optional suction.

C: Blunt trocar for removal of sample

D: Sealing cap for inner cannula (optional use)



Product #	Description	Standard Wire Gauge
72-230053	Bergstrom Needle 3mm Ø x 100mm	10 swg x 4"
72-230054	Bergstrom Needle 4mm Ø x 100mm	8 swg x 4"
72-230055	Bergstrom Needle 4.5mm Ø x 100mm	7 swg x 4"
72-230056	Bergstrom Needle 5mm Ø x 100mm	6 swg x 4"
72-230057	Bergstrom Needle 6mm Ø x 100mm	5 swg x 4"

Ask your product consultant for coordinating Sealing Caps!