



**MILLENNIUM  
SURGICAL**

**SURGICALINSTRUMENTS.COM**

UNMATCHED SERVICE • UNLIMITED INSTRUMENTS

822 Montgomery Ave #205  
Narberth, PA 19072  
(800) 600-0428

# Phlebectomy Instruments



5-47    5-710    5-716    5-5284    5-5306    5-18SC1474    5-87TC



5-80TC    5-6121    5-6125    5-1176    5-3200200    5-3800

Item Number	Description	Quantity
5-47	Knife handle no. 3 fitting surgical blades Nos. 10 through 15	1
5-710	Halstead misquito forceps 5" (12.7 CM), curved, extra delicate	6
5-716	Halstead misquito forceps 5" (12.7 CM), curved, 1x2 teeth	1
5-5284	Metzenbaum scissors 5" (12.7 CM), curved, delicate pattern	1
5-5306	Iris Scissors 4.5" (11.4 CM), curved	1
5-18SC1474	Stevens Tenotomy Scissors 4.5" (11.4 CM), curved, sharp/sharp, SuperCut	1
5-87TC	Webster Needle Holder 4.75" (12.1 CM), smooth jaws Tungsten Carbide	1
5-80TC	Crile-Wood Needle Holder 6" (15.2 CM), serrated jaws, Tungsten Carbide	1
5-6121	Adson Tissue Forceps 4.75" (12.1 CM), 1x2 teeth and cross serrated platform	1
5-6125	DeBakey-Adson Tissue Forceps 4.75" (12.1 CM), atraumatic tips 1.5 MM wide	1
5-1176	Senn Retractor 6 3/8" (16.2 CM), blunt, double ended	1
5-3200200	Tray with Silicone Mats 7" x 11" (17.8 CM x 27.9 CM), double layer	1
5-3800	Instrument care kit system spray lube (5-3700), surgical instrument cleaner (5-3720), surgical instrument stain remover (5-3740), instrument care DVD (5-3810), nylon brush and stainless steel brush	1



**MILLENNIUM  
SURGICAL**

**SURGICALINSTRUMENTS.COM**

UNMATCHED SERVICE • UNLIMITED INSTRUMENTS

822 Montgomery Ave #205  
Narberth, PA 19072  
(800) 600-0428

# Crochet Phlebectomy Hooks



Tip enlarged to show detail

Crochet hooks are popular with clinicians due to the variety of sizes available. These hooks are a solid choice for exteriorizing the vein through the incision. Crochet patterns have a blunt, rounded tip and provide optimal elevation for all removal techniques.  
(German, stainless steel)



-  5-10403 Style 00 3.00MM 6" (15.2 CM) Length
-  5-10404 Style 0 2.75 MM6" (15.2 CM) Length
-  5-10405 Style 1 2.60 MM6" (15.2 CM) Length
-  5-10406 Style 2 2.45 MM6" (15.2 CM) Length
-  5-10407 Style 3 2.30 MM6" (15.2 CM) Length
-  5-10408 Style 4 2.15 MM6" (15.2 CM) Length
-  5-10409 Style 5 2.00 MM6" (15.2 CM) Length
-  5-10410 Style 6 1.85 MM6" (15.2 CM) Length
-  5-10411 Style 7 1.70 MM6" (15.2 CM) Length
-  5-10412 Style 8 1.55 MM6" (15.2 CM) Length

5-1040 Set. Includes all 10 Styles #00-8

## Mueller-Style Hooks

-  5-10321 Right #1 4.75" (12.1 CM) Length
-  5-10322 Right #2 5" (12.7 CM) Length
-  5-10323 Right #3 5.25" (13.3 CM) Length
-  5-10324 Right #4 5.25" (13.3 CM) Length
- 5-1032 Set. Includes all 4 right hooks, sizes #1-4
-  5-10331 Left #1 4.75" (12.1 CM) Length
-  5-10332 Left #2 5" (12.7 CM) Length
-  5-10333 Left #3 5.25" (13.3 CM) Length
-  5-10334 Left #4 5.25" (13.3 CM) Length
- 5-1033 Set. Includes all 4 left hooks, sizes #1-4





**MILLENNIUM  
SURGICAL**

**SURGICALINSTRUMENTS.COM**

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205  
Narberth, PA 19072  
(800) 600-0428

# Ramlet Phlebectomy Hooks



Fine, sharp, angular hook and short, round handle. Facilitates dissection and provides a secure grip on the vein as well as moderate traction during vertical avulsion of the vein. Minimizes damage to surrounding tissues and lymphatics. Handle design allows for the instrument to be rolled between the fingers during dissection, lessening hand and wrist movements. Available in right and left handed patterns. (German, stainless steel)



5-10351	Right #1	4" (10.2 CM) Length
5-10352	Right #2	4" (10.2 CM) Length
5-1035	Set. Includes Right Hooks Sizes #1-2	
5-10361	Left #1	4" (10.2 CM) Length
5-10362	Left #2	4" (10.2 CM) Length
5-1036	Set. Includes Left Hooks Sizes #1-2	

## Oesch-Style Hooks



5-10381	Left #1	6.5" (16.5 CM) Length
5-10382	Left #2	6.5" (16.5 CM) Length
5-10383	Left #3	6.5" (16.5 CM) Length
5-10371	Right #1	6.5" (16.5 CM) Length
5-10372	Right #2	6.5" (16.5 CM) Length
5-10373	Right #3	6.5" (16.5 CM) Length

5-1037: Set, includes all 3 Right Hook sizes #1-3  
5-1038: Set, includes all 3 Left Hook sizes #1-3



**MILLENNIUM  
SURGICAL**

**SURGICALINSTRUMENTS.COM**

UNMATCHED SERVICE • UNLIMITED INSTRUMENTS

822 Montgomery Ave #205  
Narberth, PA 19072  
(800) 600-0428

# Varady Extractors

Double ended with a variety of spatula and hook designs to meet all clinical needs. Spatula options help facilitate dissection and mobilization of the vein prior to elevation for extraction. Various hook designs and sizes are available to provide controlled elevation of the varix for exteriorization prior to proximal and distal avulsion. Specialization of Varady extractors provide optimal cosmetic and clinical results. (German, stainless steel)



5-10305



5-10306



5-10307



5-10308



5-10309



5-10310

Item Number	Description	Tips	Length
5-10305	Dissecting Spatula	1.9 MM and 2.9 MM width working ends	7" (17.8 CM)
5-10306	Large Hook w/ Ball Tip & Spatula	2.8 MM wide hook and 1.9 MM wide spatula	7" (17.8 CM)
5-10307	Medium Crochet Style Hook & Spatula	1.8 MM wide hook and 1.7 MM side spatula	7" (17.8 CM)
5-10308	Micro Spatula	0.7 MM and 1.2 MM width working ends	6.75" (17.1 CM)
5-10309	Micro Hook & Spatula	1.3 MM wide hook and 0.7 MM wide spatula	6.75" (17.1 CM)
5-10310	Small Hook & Spatula	1.5 MM wide hook and 1.1 MM wide spatula	6.75" (17.1 CM)
5-1030	Set, includes all 6 Extractors		