



**MILLENNIUM
SURGICAL**


SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Bent Hohmann Retractors

The Bent Hohmann Retractor is useful for retracting tissues at the margins of the joint. In hip arthroplasty it is especially useful for retracting tissues in the posterior acetabular brim. In knee arthroplasty it can be passed over the margins of the joint and held in place with weights or by hand.



Ordering Information

Wide

Item Number 101-6590
Overall length: 9.375"
Handle length: 7"
Blade width: 41mm
Depth from bend: 4.75"

Wide w/ Extra Long Handle

Item Number 101-659001
Overall length: 11"
Handle length: 9"
Blade width: 41mm
Depth from bend: 5.5"

Narrow

Item Number 101-7110
Overall length: 9.75"
Handle length: 7"
Blade width: 19mm
Depth from bend: 4.75"

Narrow w/ Extra Long Handle

Item Number 101-711001
Overall length: 11.5"
Handle length: 10"
Blade width: 19mm
Depth from bend: 4.75"

Short-Tipped Narrow

Item Number 101-71115
Overall length: 9.25"
Handle length: 7"
Blade width: 19mm
Depth from bend: 4.25"

Short-Tipped Narrow w/ Extra Long Handle

Item Number 101-711501
Overall length: 11"
Handle length: 10"
Blade width: 19mm
Depth from bend: 4.25"