



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Large Bone Curettes

Designed with a 5/16" diameter shaft allowing better visualization into the medullary canal during use. The contoured handle is designed to keep the curette from slipping in the surgeon's hand and for better control. Shaft length: 10-1/2".

The angled large curette is designed for use in the acetabulum or exposed bone. The shaft diameter is 5/16" and has a contoured handle. Shaft length 10-1/2".



Storage Case Available



Quotations Ends at Actual Size

Ordering Information

Set of 1 each, including case	Item Number 101-516000
Angled Small	Item Number 101-516001
Straight Small	Item Number 101-516002
Angled Medium	Item Number 101-516003
Angled Large	Item Number 101-516004
Straight Medium	Item Number 101-516005
Storage Case Available!	