



MILLENNIUM
SURGICAL

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE - UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Reusable Orthopedic & Podiatric Drill Bits

Twist Drill Bits

Fits any Jacobs chuck

Drill bits with Jacob chuck end.

OAL = Overall Length

FL - Flute Length

	Diameter	OAL/FL
85-1517127	1.5mm	127-17mm
85-2023127	2.0mm	127-23mm
85-2423127	2.4mm	127-23mm
85-2523127	2.5mm	127/23mm
85-2730127	2.7mm	127/30mm
85-3242127	3.2mm	127/42mm
85-3542127	3.5mm	127/42mm
85-4034127	4.0mm	127/34mm
85-4534127	4.5mm	127/34mm
85-4734127	4.7mm	127/34mm
85-5042127	5.0mm	127/42mm



Twist Drill Rack: 98.8127

Holds 13 127mm twist drills

(drill bits sold separately)



Holds:

2 - 1.5mm

2 - 2.0mm

2 - 2.5mm

1 each: 2.7, 3.2,

3.5, 4.0, 4.5,

4.7, and 5.0mm



Twist Drill Rack: 98.8178

SQC 7" for Drill Bits SQC

w/ OD 1.10 to 3.50mm

6 3/4" x 4 1/4" x 3/4"

(drill bits sold separately)



MILLENNIUM
SURGICAL

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE - UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Reusable Orthopedic & Podiatric Drill Bits



SQC Drill Bits

	Diameter	OAL/FL
86.8211	1.1mm	60/20mm
86.8215	1.5mm	85/23mm
86.8216	1.5mm	110/23mm
86.8220	2.0mm	100/23mm
86.8221	2.0mm	125/23mm
86.8222	2.2mm	110/35mm
86.8224	2.5mm	110/30mm
86.8226	2.5mm	110/35mm
86.8225	2.5mm	180/35mm
86.8227	2.7mm	110/35mm
86.8228	2.7mm	125/35mm
86.8232	3.2mm	145/42mm
86.8233	3.2mm	195/42mm
86.8235	3.5mm	110/42mm
86.8236	3.5mm	195/42mm
86.8240	4.0mm	195/42mm
86.8245	4.5mm	145/42mm
86.8246	4.5mm	195/42mm

tin coated

SQC = Small quick coupling

OAL = Overall length

FL = Flute length

- Drill bits with quick coupling ends
- Designed to fit quick coupling handles 86.0040 and 86.0045 and power adaptor 86.1002

SQC Three-Fluted Drill Bits

	Diameter	OAL/FL
86.8725	2.5mm	230/30mm
86.8732	3.2mm	230/30mm
86.8765	4.5mm	195/45mm

SQC = Small quick coupling

OAL = Overall length

FL = Flute length

- Drill bits with quick coupling ends
- Designed to fit quick coupling handles 86.0040 and 86.0045 and power adaptor 86.1002
- 86.8725 and 86.8732 are calibrated
- 86.8765 is not calibrated

86.8725 86.8765
86.8732

SQC Four-Fluted Twist Drill Bits

	Diameter	Max ID	FL
86.8620	2.0mm	1.2mm	25mm
86.8627	2.7mm	1.3mm	25mm
86.8635	3.5mm	1.8mm	35mm
86.8640	4.0mm	2.0mm	35mm
86.8645	4.5mm	2.1mm	40mm

- 5"
- Straight Shank
- 4 Flutes
- Round end
- *Cannulated!*