



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Unger Anterior Total Hip Instruments

Universal instrument system specifically designed for
Direct Anterior approach THR

Designed by Anthony Unger, MD



Unger Anterior Total Hip Instruments

Universal instrument system specifically designed for
Direct Anterior approach THR



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Ordering and Product Information

Numbers correspond to the instruments featured on the first page. Please reference the following product number(s) when speaking with a Product Consultant.

	Product Number	Name	Description
1	101-3001	Unger Wide Hohmann - Single Prong	Blade Width: 43.6mm Blade Depth: 5" Overall Length: 13.5"
2	101-3008	Unger Wide Hohmann - Double Prong	Blade Width: 33.6mm Blade Depth: 5" Overall Length: 13.5"
3	101-3002	Unger Narrow Hohmann	Blade Width: 24.5mm Blade Depth: 4" Overall Length: 13"
4	101-3003	Unger Blunt Narrow Cobra	Blade Width: 22.5mm Blade Width at Tip: 12mm Blade Depth: 5.25" Overall Length: 14.5"
5	101-3004	Unger Canal Finder Rasp - Straight	Overall Length: 11" Handle Length: 5"
6	101-300401	Unger Canal Finder Rasp - Coved	Overall Length: 11" Handle Length: 5"
7	101-3005R	Unger Box Osteotome - Right	Overall Length: 12"
8	101-3005L	Unger Box Osteotome - Left	Overall Length: 12"
9	101-3006	Unger Femoral Neck Elevator	Blade Width at Widest: 25mm Overall Length: 12.5" Handle Length: 8.25"
10	101-3007	Unger Soft Tissue Protector	Blade Width: 50.5mm Blade Depth: 1.75" Overall Length: 10.125"
11	101-3009R	Unger Offset Narrow Hohmann - Right	Blade Width: 24.5mm Blade Depth: 4" Overall Length: 13"
12	101-3009L	Unger Offset Narrow Hohmann - Left	Blade Width: 24.5mm Blade Depth: 4" Overall Length: 13"