



**MILLENNIUM  
SURGICAL**

[SURGICALINSTRUMENTS.COM](http://SURGICALINSTRUMENTS.COM)

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205  
Narberth, PA 19072  
(800) 600-0428

# Universal Screwdriver Set

Helps eliminate the opening of multiple sterile packs when a specific size or style of screwdriver is needed

Helpful during revision total joint surgery where screws have been used, during the removal of bone plates, fracture fixation screws or bone graft screws. Set consists of: 6 (six) double ended screw driver bits — small & large single slot, cross & cruciate, 3.5mm & 4.5mm hex, small & large phillips, small & large torx — a handle which accommodates any of the above bits, and a sterilization case.



Set in Storage Case

Product Numbers	
101-5195	Complete set with storage case (also sold individually)
101-519501	Handle
101-519502	Straight (single slot) Large: 7x1.5mm. Small: 5x1mm
101-519503	Cross / Cruciate Large: 7mm. Small: 6mm
101-519504	Hex Large: 4.5mm. Small: 3.5mm
101-519505	Phillips Large: 4mm. Small: 3.5mm
101-519506	Small Torx: #10 & #15
101-519507	Large Torx: #20 & #25