



**MILLENNIUM  
SURGICAL**

**SURGICALINSTRUMENTS.COM**

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205  
Narberth, PA 19072  
(800) 600-0428

# Arthroscopes & Sinusscopes



This line of rigid scopes offers superior image. Higher light fiber density ensures more light at the operative site, while the sophisticated lens systems provide optimum performance, image quality and optical efficiency.

Featuring soldered Sapphire scratch resistant windows and innovative rigid rod lens systems, ensuring the reliability you need.

Endoscopes are Guaranteed for autoclavability for 1 year.

Seamlessly compatible with Karl Storz, Conmed and Linvatec style cannula systems.

Endoscopes are designed and manufactured in the USA!

Arthroscopes / Sinusscopes	
Part Number	Product Description
47-MVA2701	2.7 mm, 0°, 175 mm working length
47-MVA2731	2.7 mm, 30°, 175 mm working length
47-MVA400	4.0 mm, 0°, 175 mm working length
47-MVA430	4.0 mm, 30°, 175 mm working length
47-MVA445	4.0 mm, 45°, 175 mm working length
47-MVA470	4.0 mm, 70°, 175 mm working length

Feature	Benefit
Soldered eyepiece window	Eliminates leaks, provides a longer working life for the scope
Thermal expansion compensator	Assures repeatable autoclavability
Sapphire objective window and eyepiece	Eliminates scratching
Maximum fiber packaging capacity	Provides greater light at operative site - More fibers in scope
Color-coded EZ-ID™ eyepieces	Provides easy identification - No need for ID tape
US manufacturer to FDA quality system	Guarantees quality