



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Laparoscopes



Performance

Our rigid scopes offer superior imaging. Higher light fiber density ensures more light at the operative site, while our sophisticated lens systems provide optimum performance, image quality and optical efficiency.

Bring out the Best of your HD Camera

Paired with your HD Camera systems, our scopes provide superior visualization in the OR.

Durability

Featuring soldered Sapphire scratch resistant windows and innovative rigid rod lens systems, ensuring the reliability you need.

Autoclavability

MediVision endoscopes are Guaranteed for autoclavability for 1 year.

Compatibility

Seamlessly compatible with all types of laparoscopic cannulas.

Reliability

Designed and manufactured in the USA..

Product Ordering Information

Laparoscopes

MV5500	5.5 mm, 0°, 300 mm working length
MV5530	5.5 mm, 30°, 300 mm working length
MV5545	5.5 mm, 45°, 300 mm working length
MV1000	10.0 mm, 0°, 330 mm working length
MV1030	10.0 mm, 30°, 330 mm working length
MV1045	10.0 mm, 45°, 330 mm working length

Feature:	Benefit:
Soldered eyepiece window	Eliminates leaks, providing a longer working life for the scope
Thermal expansion compensator	Assures repeatable autoclavability
Sapphire objective and eyepiece window	Eliminates scratching
Maximum fiber packing capacity	Provides greater light at operative site - More fibers in scope
Color coded EZ-ID™ eyepieces	Provides easy identification - No need for ID tape
US manufacturer to FDA quality system	Guarantees quality

Service as well as repair exchange programs available for all scopes. We also offer a high quality, economical and rapid solution to your endoscopic service needs. Trademarked items are the property of their respective owners.