

Surgical Instrument Care Catalog

Brushes




**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Suction Tube Cleaning Brushes (Baron & Frazier): Pack of 3



Item# 10-1394	1 mm	6"	Pack/3
Item# 10-1396	1.6 mm	6"	Pack/3
Item# 10-1400	2 mm	8"	Pack/3
Item# 10-1401	2.3 mm	8"	Pack/3
Item# 10-1402	2.6 mm	8"	Pack/3
Item# 10-1403	3 mm	8"	Pack/3
Item# 10-1404	3.3 mm	8"	Pack/3
Item# 10-1405	3.6 mm	8"	Pack/3

Double End Valve Brushes: Pack of 3



Item# 10-1398	5.75" length 9 mm and 5 mm diameter tips Pack/3
---------------	---

Brush Drying Rack & Brush Holder



Item# 10-1399	1 Brush Rack
Item# 10-1545	1 Brush Holder



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM


UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Surgical Instrument Care Catalog

Brushes

General Instrument Cleaning Brush: Packs of 3, 12 and 24



Soft Stainless Steel		Nylon
Item# 10-1650	Pack/3	Item#: 10-1651
Item# 10-1652	Pack/3	Item#: 10-1657
Item# 10-1658	Pack/3	Item#: 10-1659

Disposable Instrument Cleaning Brush: Packs of 12 and 24



Disposable	
Item# 10-1465	Pack/12
Item# 10-1648	Pack/24

Double-End Instrument Cleaning Brush: Pack of 3



Nylon	
Item# 10-1444	Pack/13



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Surgical Instrument Care Catalog

Brushes

Stainless Steelware Cleaning Brush: Pack of 3



Use for bowls, basins and other stainless steelware.

2.5" wide x 2/5" long brush head. 5" plastic handle
White medical grade nylon bristles are 1.5" long.

Item# 10-1460 Pack/3

Medical Surface Cleaning Brush: Pack of 3



Use on sinks, autoclaves, ultrasonic cleaners,
bowls, basins, etc.

3" wide by 3.25" long brush head. 4.5" plastic handle,
white medical grade nylon bristles are 2" long.

- | | |
|--------|----------------|
| Pack/3 | Item#: 10-1418 |
| | Item#: 10-1419 |
| | Item#: 10-1421 |

Color coded brushes may be used in "critical control points" within the hospital, in a department or by utility as a means of controlling cross-contamination or simply to identify ownership.

The product has been modeled on the FDA's HACCP program for prevention of cross contamination in food borne illness.

Large Instrument Cleaning Brush: Pack of 3



Use on retractors, specula, etc.

1" wide x 5" long brush head. 5" plastic handle.
White medical grade nylon bristles are 1" long.

Item# 10-1420 Pack/3

Autoclave Cleaning Brush: Pack of 3



Use on autoclave chambers and doors, sinks, carts, etc.

3" wide by 3.5" long brush head. 16" plastic handle.
White medical grade nylon bristles are 2" long.

Item# 10-1422 Pack/3

Autoclave Cleaning Brush Kit: 1 Each



- (1) Telescopic handle extends from 24" to 48"
- (1) Large cleaning brush, 5.5" wide by 7.5" long with round corner bristles.
- (1) Hand brush, 4.5" handle.

Item# 10-1424 1 Each



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428


Instrument Stringers / Trays

E-Z Long Instrument



Item# 10-1806 6" x 2.5"
Item# 10-1808 8" x 2.5"
Item# 10-1810 10" x 2.5"
Item# 10-1812 12" x 2.5"
Item# 10-1816 16" x 2.5"

Weinstein Instrument Stringer




Item# 10-1800 5.5"
Item# 10-1805 11.2"

Bunt Instrument Stringer



Item# 10-1770

Saunders Instrument Rack



Item# 10-1780

Instrument Cath Tray



Item# 10-1738
8-7/8" x 5" x 2"

Lids sold separately.
Item# 10-1722

Instrument Tray



Item# 10-1877
12.25" x 7.75" x 2.25"

Lids sold separately.
Item# 10-1878

Sterilization Trays



Item# 10-1705
6-3/8" x 3-3/8" x 1-7/5"
NOT SHOWN

Item# 10-1710
10" x 6.5" x 2.5"

Item# 10-1720
20" x 10.5" x 3.5"

Mayo Tray



Item# 10-1691
13-3/4" x 9-3/4" x 3/4"

Item# 10-1694
19-1/8" x 12-5/8" x 3/4"



**MILLENNIUM
SURGICAL**

SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205
Narberth, PA 19072
(800) 600-0428

Non Latex Instrument Protectors

Sklartip™ Vented Instrument Protectors: Pack of 100



Item# 96-1400
11.0mm - 1/32" x 3/4"
Maroon Pack/100



Item# 96-1402
6.0mm - 3/16" x 1"
Yellow Pack/100



Item# 96-1404
4.0mm - 1/8" x 1"
Red Pack/100



Item# 96-1406
3.0mm - 7/64" x 3/4"
Green Pack/100



Item# 96-1408
2.0mm - 5/64" x 3/4"
Blue Pack/100



Item# 96-1410
1.5mm - 1/16" x 3/4"
White Pack/100



Item# 96-1412
25 x 3.0mm
1/8" x 1" x 1"
Black Pack/100



Item# 96-1414
18 x 3.0mm
5/64" x 5/8" x 1"
Brown Pack/100



Item# 96-1416
10 x 3.0mm
1/16" x 3/8" x 1"
Orange Pack/100

Sklartip™ Non-Vented Instrument Protectors: Pack of 100



Item# 96-1500
11.0mm - 1/32" x 3/4"
Maroon Pack/100



Item# 96-1502
6.0mm - 3/16" x 1"
Yellow Pack/100



Item# 96-1504
4.0mm - 1/8" x 1"
Red Pack/100



Item# 96-1506
3.0mm - 7/64" x 3/4"
Green Pack/100



Item# 96-1508
2.0mm - 5/64" x 3/4"
Blue Pack/100



Item# 96-1510
1.5mm - 1/16" x 3/4"
White Pack/100



Item# 96-1512
25 x 3.0mm
1/8" x 1" x 1"
Black Pack/100



Item# 96-1514
18 x 3.0mm
5/64" x 5/8" x 1"
Brown Pack/100



Item# 96-1516
10 x 3.0mm
1/16" x 3/8" x 1"
Orange Pack/100