


Endoscopic Sterilization Trays



2000-100-001
Camera case with brackets
ext size: 9" x 10" x 3"
int size: 8.47" x 8.68" x 2.55"




2000-100-021
Light cord case
ext size: 9" x 9" x 1.5"
int size: 8.57" x 8.72" x 1.18"




2000-100-011
Scope case 11" without brackets
ext size: 3" x 11" x 2.5"
int size: 2.93" x 10.93" x 2.5"



2000-100-019
Scope case 18" without brackets
ext size: 3" x 18" x 2.5"
int size: 2.93" x 17.93" x 2.46"



2000-100-005
Scope case, two level, with brackets/camera & light cord compartments
ext size: 9" x 20" x 4"
int size: 8.47" x 18.68" x 3.55"



2000-100-004
Scope case, two level, without brackets/dividers
ext size: 9" x 20" x 4"
int size: 8.47" x 18.68" x 3.55"



2000-100-006
Scope case for 2-4 scopes with brackets
ext size: 4" x 17.5" x 1.5"
int size: 3.57" x 17.19" x 1.35"



2000-100-035
Surgical case, without brackets/pin mat
ext size: 10" x 21.5" x 2.5"
int size: 9.72" x 20.65" x 2.42"



2000-100-015
Lap chole case, two level, with insert and brackets for up to 35 cm or 13.5" instruments (shown with optional pin mat in the base for additional instrument storage)
ext size: 11" x 23" x 8"
int size: 10.7" x 21.62" x 7.66"



2000-100-017
Lap chole case, two level, with insert and brackets for up to 65 cm or 26" instruments (shown with optional pin mat in the base for additional instrument storage)
ext size: 11" x 27" x 8"
int size: 10.7" x 25.62" x 7.66"


millenniumsurgical.com - 822 Montgomery Ave #205 - Narberth, PA 19072
www.surgicalinstruments.com - (P) 800.600.0428 - (F) 800.600.0429

Lap Chole Sterilization Trays




2000-100-026

Lap chole 14 instrument insert, 8" A-type container (fits in a 10" x 20" x 6" container)
ext size: 9.5" x 20" x 5.5"



2000-100-018

Lap chole 14 instrument insert special assembly (fits in a 10" x 20" x 6" container)
ext size: 9.5" x 20" x 5.5"




2000-100-015

Lap chole case, two level, with insert and brackets for up to 35 cm or 13.5" instruments (shown with optional pin mat in the base for additional instrument storage)
ext size: 11" x 23" x 8"
int size: 10.70" x 21.62" x 7.66"



2000-100-017

Lap chole case, two level, with insert and brackets for up to 65 cm or 26" instruments (shown with optional pin mat in the base for additional instrument storage)
ext size: 11" x 27" x 8"
int size: 10.7" x 25.62" x 7.66"



2000-100-079

Lap chole 14 instrument insert (fits in a 10" x 20" x 6" container)
ext size: 9.9" x 19.5" x 4.9"

millenniumsurgical.com - 822 Montgomery Ave #205 - Narberth, PA 19072
www.surgicalinstruments.com - (P) 800.600.0428 - (F) 800.600.0429

Osteotome Sterilization Trays



2000-100-027

Osteotome case for curved
and straight Lambottes
ext size: 9" x 11.5" x 3"
int size: 8.57" x 11.12" x 2.38"


millenniumsurgical.com - 822 Montgomery Ave #205 - Narberth, PA 19072
www.surgicalinstruments.com - (P) 800.600.0428 - (F) 800.600.0429

Rongeur Forcep Sterilization Trays



2000-100-010

Forcep case for up to 8
Rongeurs designed for size
4 & 5 handles with handle guide
ext size: 9" x 13" x 6"
int size: 8.59" x 12.09" x 5.84"



2000-100-025

Forcep case for up to 16
arthroscopic and or sinus
forceps interlacing 8 forceps in
one direction and eight in the
other direction
ext size: 9" x 15" x 5"
int size: 8.59" x 13.84" x 4.84"

millenniumsurgical.com - 822 Montgomery Ave #205 - Narberth, PA 19072
www.surgicalinstruments.com - (P) 800.600.0428 - (F) 800.600.0429