



**MILLENNIUM  
SURGICAL**

**SURGICALINSTRUMENTS.COM**

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205  
Narberth, PA 19072  
(800) 600-0428

# Ultrasonic Cleaners



	Model 1800	Model 2800	Model 3800	Model 5800	Model 8800
Tank Size (L-W-D)	6"x5.5"x4"	9.5"x5.5"x4"	11.5"x6"x6"	11.5"x9.5"x6"	19.5"x11.5"x6"
Overall Size (L-W-D)	10"x12"x11.5"	13.5"x12"x11.5"	16"x12"x14.5"	16"x15.5"x14.5"	24"x18"x14.5"
Tank Capacity (Liters)	0.5 gal. (1.9l)	.75 gal (2.8l)	1.5 gal (5.7l)	2.5 gal (9.5l)	5.5 gal (20.8l)
Weight	7 lbs	9 lbs	12 lbs	14 lbs	26 lbs
Frequency	40 kHz	40 kHz	40 kHz	40 kHz	40 kHz
Drain	NO	NO	YES	YES	YES

The above models are supplied with a tank cover. Accessories are also available, including perforated trays, solid trays, beaker covers, and support racks. Note, parts should not be allowed to rest on the tank bottom.



**MILLENNIUM  
SURGICAL**

**SURGICALINSTRUMENTS.COM**

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205  
Narberth, PA 19072  
(800) 600-0428

# Ultrasonic Cleaners

## Model 1800

Part Number	Model	Description
CPX-952-116R	M1800	Mechanical Timer
CPX-952-117R	M1800H	Mechanical Timer and Heater
CPX-952-118R	CPX1800H	Digital Timer, Heater and Degas, Temp Monitor
CPX-952-119R	CPX1800	Digital Timer
100-032-515		Tank Cover (included with unit)
100-410-170		Solid Tray
100-410-160		Perforated Tray
100-032-516		Beaker Cover (1-600ml)
100-916-333		Mesh Basket
CPN-916-039		Support Rack

## Model 1800

Part Number	Model	Description
CPX-952-316R	M3800	Mechanical Timer
CPX-952-317R	M3800H	Mechanical Timer and Heater
CPX-952-318R	CPX3800H	Digital Timer, Heater and Degas, Temp Monitor
CPX-952-319R	CPX3800	Digital Timer
100-032-519		Tank Cover (included with unit)
100-410-174		Solid Tray
100-410-164		Perforated Tray
100-032-520		Beaker Cover (2-600ml)
100-032-528		Beaker Cover (3-250ml)
100-916-335		Mesh Basket
CPN-916-041		Support Rack

## Model 2800

Part Number	Model	Description
CPX-952-216R	M2800	Mechanical Timer
CPX-952-217R	M2800H	Mechanical Timer and Heater
CPX-952-218R	CPX2800H	Digital Timer, Heater and Degas, Temp Monitor
CPX-952-219R	CPX2800	Digital Timer
100-032-517		Tank Cover (included with unit)
100-410-172		Solid Tray
100-410-162		Perforated Tray
100-032-518		Beaker Cover (2-600ml)
100-032-527		Beaker Cover (2-250ml)
100-916-334		Mesh Basket
CPN-916-040		Support Rack

## Model 5800

Part Number	Model	Description
CPX-952-516R	M5800	Mechanical Timer
CPX-952-517R	M5800H	Mechanical Timer and Heater
CPX-952-518R	CPX5800H	Digital Timer, Heater and Degas, Temp Monitor
CPX-952-519R	CPX5800	Digital Timer
100-032-510		Tank Cover (included with unit)
100-410-176		Solid Tray
100-410-166		Perforated Tray
100-032-511		Beaker Cover (4-600ml)
CPN-916-042		Support Rack

## Model 8800

Part Number	Model	Description
CPX-952-116R	M8800	Mechanical Timer
CPX-952-117R	M8800H	Mechanical Timer and Heater
CPX-952-118R	CPX8800H	Digital Timer, Heater and Degas, Temp Monitor
CPX-952-119R	CPX8800	Digital Timer
100-917-050		Tank Cover (included with unit)
100-410-178		Solid Tray
100-410-168		Perforated Tray
CPN-246-014		Beaker Cover (6-600ml)
CPN-916-043		Support Rack

Additional accessories available! Ask your Instrument Expert for more info!



The above models are supplied with a tank cover. Accessories are also available, including perforated trays, solid trays, beaker covers, and support racks. Note, parts should not be allowed to rest on the tank bottom.